BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFIA/BIBLIOGRAFIA

1. Henderson, A.R., Donald W.M., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 5th Ed., Butris, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 352.
2. Scherwin, J.E., *Liver function Clinical Chemistry Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A. Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
3. German Society for Clinical Chemistry, *Z. Klin. Chem. Klin. Biochem.*, (1972), **10**, 281.
4. Committee of Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology, *Scand. J. Clin. Lab. Invest.*, (1974), **33**, 291.
5. Guder, World Health Organization, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*, WHO/DIL/LAB/99.1 Rev.2, (2002), 21.
6. Rosaki, S.B., *Clin. Chem.*, (1993), **37**, 648.
7. Vassault A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
8. Vassault A., et al., *Ann. Biol. Clin.*, (1999), **57**, 685.
9. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
10. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
11. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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English - EN
INTENDED USE

ELITech Clinical Systems ALP (DEA) SL is an *in vitro* diagnostic reagent intended for the quantitative determination of alkaline phosphatase in human serum samples.

CLINICAL SIGNIFICANCE (1-2)

Alkaline phosphatase (ALP) corresponds to a group of phosphatases that display maximum activity at alkaline pH. ALP is widely distributed in liver, osteoblasts, intestinal epithelium, kidneys, and placenta.

The rate of ALP rises physiologically for children and teenagers during periods of active growth, as well as for women during the third trimester of pregnancy. Marked increases of ALP rate are observed in case of extra-hepatic obstruction (gallstones, tumours...) and bone diseases such as Paget's disease and osteogenic bone cancer.

ALP activity can also increase moderately in the event of intra-hepatic obstruction, hepatitis, cirrhosis, or with rickets, osteomalacia, hyperparathyroidism, healing of bone fractures.

M

ALP (DEA) SL

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PASL-0230	R1 4 x 20 mL	R2 4 x 5 mL
PASL-0400	R1 2 x 50 mL	R2 1 x 26 mL
PASL-0420	R1 4 x 50 mL	R2 2 x 26 mL

PRODUCT DETERIORATION
 - The product should be clear. Cloudiness would indicate deterioration.
 - Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
 - Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES (1,2,5)

Specimen
 - Serum free from hemolysis.
 - Do not use other specimens.
Warnings and precautions
 - According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.
 - It is better analyse freshly specimen (no later than 4 hours after collection) and to keep them at room temperature.
Storage and stability
 - ALP activity may increase when refrigerated (-2°C) or frozen (-20°C) specimen are put at room temperature.
 - If specimens are frozen for extended storage, move them to room temperature 18 to 24 hours before analysis to fully activate the enzyme.
 - Samples are stable 1 week at room temperature, 1 week at 2-8°C and 2 months at -20°C.

REFERENCE VALUES (6)

Serum (37 °C):
 Men < 270 U/L
 Women < 240 U/L

Reference values for children and teenagers during bone growth are higher than for adults.

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

Manual Procedure
 Wavelength : 405 nm
 Optical path : 1 cm
 Sample/reagent ratio : 1:50
 Temperature : 37 °C
 Read against distilled water.

Working reagent (4 volumes of R1 + 1 volume of R2) 1000 µL
Sample 20 µL

Mix and after 1 minute incubation, read absorbance at 1 minute intervals during 3 minutes. Calculate the change of absorbances per minute (ΔA/min)

Automatic Procedure
 These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

Activity (U/L) = ΔA /min x 2750

Conversion factor : U/L x 0.0167 = µkat/L

CALIBRATION

For calibration, multiparametric calibrator ELICAL 2 must be used. Its value is traceable to the manual measurement.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra XL
- Measuring range
 The reagent is linear from 20 to 900 U/L.

- Limit of Detection (7)
 Determined according to SFBC protocol, the detection limit is equal to 6 U/L.

Precision

Within-run reproducibility	n	Mean	CV (%)
		U/L	
Level 1	20	42	2,9
Level 2	20	146	0,8
Level 3	20	739	0,6

- Correlation
 A comparative study has been performed on analyser ELITech Clinical Systems Selectra XL between one-reagent procedure and two-reagent procedure on 30 serum samples. The values were between 17 and 886 U/L. The linear regression parameters are as follows:

Correlation coefficient: (r) = 0,9999
 Linear regression: $y = 0,9800 x + 3,19 \text{ U/L}$

Limitations/ Interferences (7,8)

- Do not report results outside of the usable range.
 - According to SFBC recommendations, studies have been performed to determine the level of interference from different compounds:

Unconjugated bilirubin: No significant interference up to 36 mg/L (615,8 µmol/L).
Conjugated bilirubin: No significant interference up to 25 mg/DL (427,6 µmol/L).
Haemoglobin: No significant interference up to 500 mg/dL (5 g/L).
Turbidity: No significant interference up to 600 mg/dL (6,78 mmol/L) Triglycerides equivalent.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽⁹⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽¹⁰⁻¹¹⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

On board stability/Calibration frequency

On Board Stability: 14 days
Calibration frequency: 7 days
 Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

These performances have been obtained using ELITech Selectra ProXL analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.

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