PIT-AMSL-4-v19 (06/2019)

Français - FR

Code technique : DL

☛USAGE PRÉVU

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English - EN

☛INTENDED USE

ELITech Clinical Systems AMYLASE SL is an *in vitro* diagnostic reagent intended for the quantitative determination of amylase in human serum and plasma samples.

CLINICAL SIGNIFICANCE ⁽¹⁻²⁾


α-Amylase is an enzyme from pancreatic or salivary origin that hydrolyses 1,4-α-glucosidic bonds, thus helping for starch digestion. Analysis of serum amylase is mainly used in the diagnosis of the pancreatic diseases (acute or chronic pancreatitis and their complications, carcinoma). During acute pancreatitis, a transitory increase of serum amylase is observed, a peak corresponding to of 4 to 6-fold elevation being reached in 12 to 72h after the beginning, the activity returning to the normal after 3 to 5 days. However, a serum amylase increase is also observed in other intra-abdominal pathologies, renal insufficiency, ovary cancers, salivary gland lesions, acute alcoholism, renal insufficiency or macroamylasemia (presence of a complex amylase-IgG not filtered by the glomerulus).

METHOD ⁽³⁾

Substrate: CNP-G₃ (2-chloro-4-nitrophenyl-α-maltotrioside)
 Enzymatic, Kinetic

PRINCIPE ⁽³⁾

Substrate CNP-G₃ is hydrolyzed by catalytic action of α-amylase to produce CNP (2-chloro-4-nitrophenol).

CNP-G₂ = 2-Chloro-4-nitrophenyl-α-maltoside

The rate of increase in absorbance is measured at 405 nm and is directly proportional to the activity α-amylase in the sample.

☛COMPOSITION

Reagent : R	
MES buffer, pH 6,15	50 mmol/L
Sodium chloride	70 mmol/L
Calcium chloride	6 mmol/L
Potassium thiocyanate	900 mmol/L
CNP-G ₃	2,27 mmol/L
Sodium azide	< 0,1 %

☛MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- General Laboratory equipment.
- Normal saline solution (NaCl 9 g/L)
- Do not use materials that are not required as indicated above.

☛WARNINGS AND PRECAUTIONS

- This *in vitro* diagnostic device is for professional use only.
- Contact with acids liberates toxic gas.
- The reagent contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Saliva and sweat contain amylase. It is therefore recommended to wear gloves and a mask to avoid the contamination of the reagent.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

☛STABILITIES

Store at 2-8 °C and protect from light. Do not freeze.
 Do not use after expiration dates indicated on the vial labels.
On board stability :
 The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent is ready to use.

☛PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES ⁽⁴⁾

- Specimen**
- Serum
- Lithium heparinized plasma
- Do not use other specimens.

Warnings and precautions
 According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.

Storage and stability
 Samples are stable 1 week at room temperature, 1 week at 2-8 °C and 1 year at - 20 °C.

REFERENCE VALUES ⁽⁵⁾

Sérum, plasma (37 °C): 31-107 U/L

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

For *ELITech Clinical Systems Selectra Analyzers*, applications are available on request.
Wavelength 405 nm
Temperature: 37 °C
Read against reagent blank.

Reagent R	300 µL
Sample	3 µL

Mix and after 50 second of incubation, measure the change of absorbance per minute (ΔA/min.) during 159 seconds.

With Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

$$\frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Calibrator}}} \times n \quad n = \text{Calibrator activity}$$

Conversion factor : U/L x 0.0167 = µkat/L

CALIBRATION (5)

For calibration, multiparametric calibrator ELICAL 2 must be used. Its value is traceable to IFCC method.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

☛ Measuring range

Determined according to CLSI⁽⁶⁾ EP6-A protocol, the measuring range is from 20 to 1500 U/L (0.33 to 25.00 µkat/L).

Samples having greater concentrations should be diluted 1:10 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 15000 U/L (250.00 µkat/L).

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

- Limit of Detection (LoD) and Limit of Quantification (LoQ)

Determined according to CLSI EP17-A protocol⁽⁷⁾.

$$\text{LoD} = 6 \text{ U/L (0.10 } \mu\text{kat/L)}$$

$$\text{LoQ} = 13 \text{ U/L (0.22 } \mu\text{kat/L)}$$

- Precision

Determined according to CLSI EP5-A2 protocol⁽⁸⁾.

	Mean			Within-run	Total
	n	U/L	µkat/L		
Low level	80	82	1.37	1.3	2.7
Medium level	80	204	3.40	0.9	2.2
High level	80	992	16.53	1.5	2.6

- Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (IFCC method) on 100 human serum samples according to CLSI⁽⁹⁾ EP9-A2 protocol. The values were between 21 and 1439 U/L (0.35 and 23.98 µkat/L).

The parameters of the linear regressions are as follows :

$$\text{Correlation coefficient: (r)} = 0.999$$

$$\text{Linear regression: } y = 0.976 \times x - 1 \text{ U/L (0.02 } \mu\text{kat/L)}$$

- Limitations/Interferences

- Do not report results outside of the usable range.
- Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol⁽¹⁰⁾.

Recovery is within ±10% of initial value of amylase activity of 80 and 1000 U/L.

Triglycérides: No significant interference up to 3000 mg/dL (33.9 mmol/L).
Unconjugated bilirubin: No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated bilirubin: No significant interference up to 29.5 mg/dL (504 µmol/L).
Hemoglobin: No significant interference up to 500 mg/dL.
Ascorbic acid: No significant interference up to 20.0 mg/dL.
Acetylsalicylic acid: No significant interference up to 200 mg/dL.
Acetaminophen : No significant interference up to 30 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽¹¹⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽¹²⁻¹³⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- On board stability/Calibration frequency

On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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