CK NAC SL

Références/References/
Referencias/ Referências:
CKSL-0230 4 x 20 mL + R1 4 x 5 mL
CKSL-0410 2 x 62.5 mL R1 2 x 50 mL + R2 1 x 26 mL
CKSL-0430 4 x 62.5 mL R1 4 x 50 mL + R2 2 x 26 mL

Composition du coffret/ Kit composition/
Composición del kit/ Conteúdo da embalagem :
R1 4 x 20 mL + R2 4 x 5 mL
R1 2 x 50 mL + R2 1 x 26 mL
R1 4 x 50 mL + R2 2 x 26 mL

BIBLIOGRAPHY

- Panteghini, M. & Bais R., Enzymes. Tietz Fundamentals of Clinical Chemistry, 6th Ed, Burris, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 317
- Sanhai, W.R. & Christenson R.F., Protein Isoforms: Isoenzymes and Isoforms. Clinical Chemistry: Theory Analysis, Correlation, 5th Ed, Kaplan L.A., Pesce, A.J., (Mosby, Inc.), (2010), 677 and appendix.
- Apple, F. S., Cardiac function, Tietz Fundamentals of Clinical Chemistry, 6th Ed, Burris, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 614.
- Boot, E., de Haan, L., Massive increase in serum creatine kinase during olanzapine and quetiapine treatment, not during treatment with clozapine. *Psychopharmacology* (Berlin) (2000);**150**:347-348.
- Schumann, G., et al. *Clin Chem Lab Med.*, (2002), **40**, 635-642.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature. *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
- Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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STABILITY OF REAGENTS
Store at 2-8 °C and protect from light. Do not freeze.
Do not use after expiration dates indicated on the vial labels.
On board stability:
The on-board stability is specific for each analyzer.
(Refer to § PERFORMANCE DATA).

PREPARATION
The reagents are ready to use.

REAGENT DETERIORATION
- The reagent solution should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

DAMAGED PACKAGING
Do not use the reagent if the damages of packaging might have an effect on the product performance (leakages, pierced vial).

SAMPLES (1,5)
Specimen
- Serum free from hemolysis (specimen recommended by IFCC).
- Heparinized plasma free from hemolysis.
- Do not use other specimens.

Warnings and precautions
According to Good Laboratory Practice, venipuncture should be performed prior to the administration of drugs.

Storage and stability
Samples must be analyzed immediately or stored protected from air and light for 8 hours at room temperature, 2 days at 2-8 °C or 1 month at -20 °C.

REFERENCE VALUES (5)

Men < 171 U/L
Women < 145 U/L

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

For ELITech Clinical Systems Selectra Analyzers, applications are available on request.
Wavelength: 340 nm
Temperature: 37 °C
Read against reagent blank.

Reagent R1	240 µL
Sample	12 µL

Mix and wait 4 minutes and 43 seconds, then add:

Reagent R2	60 µL
------------	-------

Mix and after a 130 second incubation, measure the change of absorbance per minute (ΔA/min.) during 106 seconds.

With Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

$$\Delta A \text{ Sample} \times n = \text{calibrator concentration}$$

$$\Delta A \text{ Calibrador}$$

$$\text{Conversion factor : } U/L \times 0,0167 = \mu\text{kat/L}$$

CALIBRATION

For calibration, multiparametric calibrator Elical 2 must be used. Its value is traceable to IFCC method⁽⁵⁾.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I (normal control) and ELITROL II (abnormal control) should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT
Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

- Measuring range
Determined according to CLSI EP6-A protocol⁽⁶⁾, the measuring range is from 10 to 1714 U/L (0,17 to 28,57 µkat/L). Samples exceeding 1714 U/L should be diluted 1:10 with NaCl 9 g/L solution (normal saline) and reassayed. Use of this procedure extends the measuring range from 1714 to 17140 U/L (28,57 to 285,67 µkat/L).

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

- Limit of Detection (LoD) and Limit of Quantification (LoQ)
Determined according to CLSI EP17-A protocol⁽⁷⁾.

$$\text{LoD} = 1 \text{ U/L (0,02 } \mu\text{kat/L)}$$

$$\text{LoQ} = 5 \text{ U/L (0,08 } \mu\text{kat/L)}$$

- Precision
Determined according to CLSI EP5-A2 protocol⁽⁸⁾.

	Mean	Within-run	Total
n	U/L	µkat/L	CV (%)
Low level	80	147	2,45
		0,7	1,7
Medium level	80	406	6,77
		1,1	2,4
High level	80	1154	19,23
		1,1	3,9

Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-cleared system equipment (IFCC method) on 100 human serum samples according to CLSI EP9-A2 protocol⁽⁹⁾. The sample concentrations were between 11 and 1712 U/L (0,18 - 28,53 µkat/L).

The parameters of the linear regressions are as follows:
Correlation coefficient: (r) = 0,998
Linear regression: $y = 1,012x + 2 \text{ U/L (0,03 } \mu\text{kat/L)}$

Limitations/Interferences

- Hemolyzed samples should not be used since significant hemolysis may lead to falsely increased CK concentration because of adenylate kinase release.

- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol⁽¹⁰⁾. Recovery within ± 10% of initial value of CK activity of 150 and 1200 U/L.

Triglycerides: No significant interference up to 3000 mg/dL (33,9 mmol/L).

Unconjugated bilirubin: No significant interference up to 30,0 mg/dL (513 µmol/L).

Conjugated bilirubin: No significant interference up to 29,5 mg/dL (504 µmol/L).

Ascorbic acid: No significant interference up to 20,0 mg/dL.

Acetylsalicylic acid: No significant interference up to 200 mg/dL.

Acetaminophen: No significant interference up to 30 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽¹¹⁾

- Many other substances and drugs may interfere. Some of them are listed in Young.⁽¹²⁻¹³⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

On board stability/Calibration frequency

On Board Stability: 28 days

Calibration frequency: 28 days

Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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