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BIBLIOGRAPHIE/BIBLIOGRAPHY/ BIBLIOGRAFÍA/BIBLIOGRAFIA

- Johnson, A.M., *Amino Acids and Proteins. Tietz Fundamentals of Clinical Chemistry*, 6th Ed, Burtis, C.A. & Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds), (2008), 286.
- Sanhai, W.R., *et al., Cardiac and Muscle Disease. Clinical Chemistry: Theory Analysis. Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc. eds.), (2010), 677 and appendix.
- Wu, A.H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 190.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCPress, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCPress, (1995).

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English - EN

INTENDED USE

ELITech Clinical Systems CRP IP is an *in vitro* diagnostic reagent intended for the quantitative determination of C-reactive protein in human serum samples on analyzers or semi-automatic analyzers. This *in vitro* diagnostic device is for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

C-reactive protein (CRP) is a non-specific protein of the acute phase of the inflammation process. Present at very low concentrations in healthy persons, CRP concentration is increased in numerous pathologies such as infections, inflammatory diseases or trauma, myocardial infarction and tumors.

In clinical practice, CRP is indicated to help in the diagnosis and follow-up of an infection or an acute or chronic inflammation.

LIMITATION OF USE

The quantitative assay of C-reactive protein alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history. Not intended for an evaluation of cardiovascular risk (CRP < 10 mg/L).

METHOD & PRINCIPLE

Immuno-turbidimetry - End Point.

The formation of CRP / anti-CRP antibody complexes is started by the addition of the antiserum to the sample in presence of an accelerator. These complexes agglutinate leading to an increase of turbidity measured at 340 nm.

COMPOSITION

Reagent 1: R1
Buffer, pH 7.43
Accelerator
Sodium azide < 0.1% (w/w)
Reagent 2: R2
Buffer, pH 7.43
Polyclonal anti-human CRP antibody (goat)
Sodium azide < 0.1% (w/w)

MATERIALS REQUIRED BUT NOT PROVIDED

- ICRP-0043 CRP IP CALIBRATOR SET
- ICRP-0046 CRP IP CONTROL I
- ICRP-0047 CRP IP CONTROL II
- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- The reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide build-up.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.

- Consult Safety Data Sheet (SDS) for a proper handling.

STABILITY
Store at 2-8 °C and protect from light. Do not freeze.
Do not use after expiration dates indicated on the vial labels.

On board stability :
The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES Specimen ⁽²⁾
- Serum
- Using any other specimen type should be validated by the laboratory.
Warnings and precautions
Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

Storage and stability ⁽⁴⁾
- 11 days at room temperature.
- 2 months at 2-8°C.
- 3 years at -20°C.

REFERENCE VALUES ⁽²⁻³⁾			
<i>Serum</i>	mg/L		mg/dL
Adults	≤ 10		≤ 1.0

Note : *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

PROCEDURE			
<u>Manual Procedure</u>			
<i>Wavelength</i> :	<i>340 nm</i>		
<i>Optical path</i> :	<i>1 cm</i>		
<i>Sample/ Reagent ratio</i> :	<i>1:17</i>		
<i>Temperature:</i>	<i>37 °C</i>		
	CALIBRATION	TEST	
Reagent R1	1000 µL	1000 µL	
Calibrator	64 µL	-	
Sample	-	64 µL	

Mix and read the absorbances (A) after an incubation of 5 minutes.(sample blank) then add :

Reagent R2	100 µL
------------	--------

Mix and read the absorbance (A2) after 4 minutes incubation.

Automatic Procedure
These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request.

CALCULATION

CRP concentration is calculated from a calibration curve obtained from five calibrators of different levels and a zero point.
Conc = f (ΔA) = f (A2 - A1)

Conversion factor : mg/L x 0.1 = mg/dL

CALIBRATION
Calibrators from CRP IP CALIBRATOR SET are traceable to ERM-DA474/IFCC reference material.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL
It is recommended that quality control sera such as CRP IP CONTROL I and CRP IP CONTROL II be used to monitor the performance of the assay. Controls have to be performed :
- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES

These performances have been obtained on Selectra Pro analyzers or earlier Selectra analyzers.

- **Measuring range**
5.0 - 190.0 mg/L (0.50 - 19.00 mg/dL)
Do not report results outside the measuring range.

Limit of Detection (LoD)

LoD =1.5 mg/L(0.15 mg/dL)

- **Hook effect**
No hook effect up to 800.0 mg/L (80.00 mg/dL).

- Precision

Representative results are presented in the following table.

		Mean		CV (%)
	n	mg/L	mg/dL	
Low level	20	16.8	1.68	2.4
Medium level	19	53.4	5.34	3.8
High level	20	97.1	9.71	1.9

		Mean		CV (%)
	n	mg/L	mg/dL	
Low level	20	19.3	1.93	5.4
Medium level	20	37.3	3.73	3.8
High level	20	105.0	10.5	1.9

- Correlation

A comparative study has been performed between CRP IP reagent on a Selectra analyzer and a similar commercially available system on 29 human serum samples.

The sample concentrations ranged from 2.6 to 211.8 mg/L (0.26 - 21.18 mg/dL). The results are as follows :
Correlation coefficient : (r) = 0.9974
Linear regression:
y = 0.9060x + 3.3 mg/L (0.33 mg/dL).

- Limitations/Interferences

- Studies have been performed to determine the level of interference from different compounds.
The following C-reactive protein levels were tested: 6 mg/L and 40 mg/L. No significant interference is defined by a recovery ±2.4 mg/L of the initial value at CRP concentration of 6 mg/L and ±15 % of initial value at CRP concentration of 40 mg/L.
Unconjugated bilirubin : No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated bilirubin : No significant inter-

ference up to 29.5 mg/dL (505 µmol/L).
Turbidity : No significant interference up to 250 mg/dL triglyceride equivalent (2.8 mmol/L).
Hemoglobin : No significant interference up to 500 mg/dL.

- Do not use turbid samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽⁵⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽⁶⁻⁷⁾

- **On board stability/Calibration frequency**
On Board Stability: 21 days
Calibration frequency: 21 days
Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

These performances have been obtained on Selectra Pro analyzers or earlier Selectra analyzers. *Results may vary if a different instrument or a manual procedure is used.*
The performances of applications not validated by ELITech are not warranted and must be defined by the user.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

TECHNICAL ASSISTANCE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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Contacte a su distribuidor local o con ELITech Clinical Systems SAS (CCsupport@elitechgroup.com)