PIT-CRCO-4-v24 (02/2019)

Français - FR

Code technique : RL

☞USAGE PRÉVU

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English - EN

☞INTENDED USE

ELITech Clinical Systems CREATININE JAFFE is an *in vitro* diagnostic reagent intended for the quantitative determination of creatinine in human serum samples.

CLINICAL SIGNIFICANCE ⁽¹⁻²⁾

Creatinine is the waste spontaneous product of creatine metabolism. It is an excellent marker of the renal function. The serum creatinine rate tends to remain constant. A high serum creatinine rate (associated to a high urea rate) corresponds to a decrease in renal glomerular filtration (FGR). The serum creatinine test is more reliable than the urea test. Indeed, the urea serum rate is affected by factors such as diet, dehydration degree and protein metabolism (the serum creatinine rate is not influenced by these factors). The test of creatinine clearance can also be used to measure the FGR. In the case of renal transplantation, any increase in serum creatinine, as little as it may be, can reflect the rejection of the transplant.

An increase of creatinine serum can be the sign of muscular necrosis.

METHOD

Colorimetric, Jaffe- Kinetic.

☞PRINCIPLE ⁽³⁻⁴⁾

The rate of formation of a coloured complex between creatinine and alkaline picrate is measured. The complex formation rate is proportional to the concentration of creatinine in the sample. The effect of interfering substances are reduced using the kinetic procedure.

☞COMPOSITION

Reagent 1: R1 Picric acid	8.73 mmol/L
Reagent 2 : R2 Sodium hydroxide	312.5 mmol/L
Disodium phosphate	12.5 mmol/L
Standard : Std Creatinine	2 mg/dL 177 µmol/L

☞MATERIALS REQUIRED BUT NOT PROVIDED

- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

☞WARNINGS AND PRECAUTIONS

- These *in vitro* diagnostic devices (reagent and standard) are for professional use only.
- The reagent R2 is classified as hazardous :

WARNING : May be corrosive to metals. Causes skin irritation. Causes serious eye irritation. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection.

- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- If eye irritation persists: Get medical advice/attention.
- Absorb spillage to prevent material damage.
- For more information, refer to the Safety Data Sheet (SDS).
- Take normal precautions and adhere to good laboratory practice.
- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.

☞STABILITIES

Store at 2-25 °C and protect from light. Do not freeze.

Do not use after expiration dates indicated on the vial labels.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents and standard are ready to use.

☞PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES ⁽⁵⁾

Specimen
- Serum

- Do not use other specimens.

☞Warnings and precautions
According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.

Storage and stability
Serum is stable for 7 days at 15 - 25 °C, 7 days at 2-8 °C and 3 months at - 20 °C.

REFERENCE VALUES ⁽⁶⁾

	Men	Women
Serum :	0.8 - 1.3 71 - 115	0.6 - 1.2 mg/dL 53 - 106 µmol/L

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

Manual Procedure	
<i>Wavelength :</i>	505 nm
<i>Optical path :</i>	1 cm
<i>Sample/reagent ratio :</i>	1:10
<i>Temperature :</i>	37 °C
<i>Read against distilled water.</i>	
Working reagent (1 volume of R1 + 1 volume of R2)	1000 µL
Sample or Standard	100 µL

Mix and read the absorbance (A1) 10 seconds after the sample or standard addition. Exactly 2 minutes after the first reading, take second reading (A2).

Automatic Procedure

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the bundle available at the end of this insert.

CALCULATION

(A2 - A1) Sample x n n = standard concentration
(A2 - A1) Standard

Conversion factor : mg/dL x 88.40 = µmol/L
mg/dL x 10 = mg/L

CALIBRATION

Concentration value of Creatinine Standard 2 mg/dL is traceable to the reference method ID-MS (Isotope Dilution - Mass Spectrometry)

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra E Analyzers

- Measuring range

The reagent is linear from 0.5 to 15 mg/dL (44 to 1326 µmol/L).

- Detection limit ⁽⁷⁾

Determined according to SFBC protocol, the detection limit is equal to 0.1 mg/dL (9 µmol/L).

- Precision

	Within-run reproducibility				Between-run reproducibility			
	n	Mean		CV %	n	Mean		CV %
		mg/dL	µmol/L			mg/dL	µmol/L	
Level 1	20	0.57	50	2.7	20	0.59	52	4.1
Level 2	20	1.60	141	0.4	20	1.60	141	2.8
Level 3	20	6.57	581	1.0	20	6.72	594	2.3

- Correlation

A comparative study has been performed on this reagent between Elitech Clinical Systems Selectra E and Cobas Mira analysers (Jaffe method) on 30 human serum samples.

The values were between 0.47 and 18.26 mg/dL (42 and 1614 µmol/L).

The parameters of linear regression are as follows :

Correlation coefficient : (r) = 0.9998

Linear regression : y = 0.9786 + 0.06 mg/dL (5 µmol/L)

- Limitations/Interferences ⁽⁷⁾

- Do not report results outside of the usable range.

- According to SFBC recommendations, some studies have been performed to determine the level of interference from different components :

Conjugated bilirubin: Negative bias from 40 µmol/L (2.3 mg/L).

Unconjugated bilirubin: Negative bias from 50 µmol/L (2.9 mg/dL).

Haemoglobin : Negative bias from 150 mg/dL.

Turbidity : No significant interference up to 580 mg/dL (6.55 mmol/L) Triglycerides equivalent .

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽⁸⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽⁹⁻¹⁰⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- On board stability/Calibration frequency

On-board stability : 3 days

Calibration frequency : 3 days

Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

These performances have been obtained using ELITech Selectra E analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.

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