PIT-IHAP-3-v8 (10/2019)

## Français - FR

### ❖ USAGE PRÉVU

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## English - EN

### ❖ INTENDED USE

ELITech Clinical Systems HAPTOGLOBIN IP is an *in vitro* diagnostic reagent intended for the quantitative determination of haptoglobin (Hp) in human samples of serum.

### CLINICAL SIGNIFICANCE (1-3)

Haptoglobin (Hp) is an acute-phase protein, which concentration increases in inflammatory syndromes.

Haptoglobin is also a protein of transport. During an intravascular hemolysis, free haemoglobin binds haptoglobin to form Hp-Hb complexes which are rapidly removed from the circulation. As a consequence Hp level is decreased. However, a decrease due to hemolysis may be compensated by an increase of Hp rate in case of concomitant inflammation. Hp level is usually decreased in hepatic diseases.

Haptoglobin dosage may be performed to appreciate inflammatory processes, and a possible hemolysis *in vivo*. Haptoglobin levels should be interpreted in conjunction with orosomucoid.

### METHOD

Immuno-turbidimetry.  
End-Point.

### PRINCIPLE

The formation of haptoglobin / anti-haptoglobin antibody complexes, started by the addition of the antiserum to the sample, is accelerated in presence of polyethylene glycol. These complexes precipitate leading to an increase of turbidity measured at 340nm. The haptoglobin concentration is determined by means of a non linear calibration curve.

### ❖ COMPOSITION

<b>Reagent 1 : R1</b>	
Phosphate Buffer	pH 7.43
Polyethylene glycol	60 g/L
Sodium azide	< 0.1 %
<b>Reagent 2 : R2</b>	
Phosphate buffer	pH 7.43
Polyclonal anti-human haptoglobin antibody (goat)	
Sodium azide	< 0.1 %

### MATERIALS REQUIRED BUT NOT PROVIDED

- Calibrator and Control :

References	Description
IPRO-0043	Protein IP Calibrator Set
CONT-0060	ELITROL I

These references are sold separately

- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

### ❖ WARNINGS AND PRECAUTIONS

- This *in vitro* diagnostic device is for professional use only.
- The reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

### ❖ STABILITIES

**Store at 2-8 °C and protect from light. Do not freeze.**  
Do not use after expiration dates indicated on the vital labels.

### PREPARATION

The reagents are ready to use.

### ❖ PRODUCT DETERIORATION

- The reagents R1 and R2 should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

- Do not use the reagent if the damages of packaging might have an effect on the product performance (leakages, pierced vial).

### SAMPLES (2,3)

**Specimen**  
- Serum free of hemolysis.  
- Do not use other specimens.  
**Warnings and precautions**  
According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.

**Storage and stability**  
Samples must be analyzed immediately or may be stored for 3 days at 2-8 °C or 2 weeks at -20°C (if frozen once).

### REFERENCE VALUES (4)

30-200 mg/dL (0,3-2,0 g/L)

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

### PROCÉDURE

This reagent can be used on most analysers, semi-automated analysers and manual method.

The applications are available on request.

Wavelength : 340 nm

Temperature : 37 °C

No sample dilution

	CALIBRATION	TEST
Reagent R1	280 µL	280 µL
Calibrator	2 µL	-
Sample	-	2 µL

Mix and read the absorbances(A1) after 75 secondes of incubation (sample blank) then add :

Reagent R2	40 µL
------------	-------

Mix and read the absorbance (A2) after 375 secondes of incubation.

### CALCULATION

Haptoglobin concentration is calculated from a calibration curve obtained from five calibrators of different levels and a zero point.

$$\text{Conc} = f(\Delta A) = f(A2 - A1)$$

### CALIBRATION

For calibration, use ELITech Protein IP Calibrator set , ref IPRO-0043, ready-to-use. Add a zero point (NaCl 9 g/L solution)

Concentration values for calibrators are traceable to ERM-DA470k/IFCC reference material.

It is recommended to make a new calibration when reagent lots change and when quality control results fall outside the established range.

### QUALITY CONTROL

To check the accuracy of assays, control serum such as ELITROL I should be used.

This control must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

### WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

### PERFORMANCE DATA at 37 °C

#### - Measuring range

The reagent is linear from 2 to 400 mg/dL (0.02 to 4.00 g/L). The exact range depends on the value of the calibrator used.

#### - Hook effect

No risk of error linked to hook effect up to 2000 mg/dL (20,00 g/L)

#### - Precision

Within-run reproducibility	n	Mean		CV (%)
		mg/dL	g/L	
Low level	20	39	0.39	1.9
Medium level	20	100	1.00	1.4
High level	20	200	2.00	2.3

Between-run reproducibility	n	Mean		CV (%)
		mg/dL	g/L	
Low level	21	34	0.34	2.4
Medium level	27	103	1.03	2.0
High level	21	143	1.43	2.0

#### - Correlation

A comparative study was performed between ELITech HAPTOGLOBIN IP and another commercial reagent (nephelometric method) on human serum samples.

The parameters of linear regression are as follows :

Correlation coefficient : (r) = 0.9945

Linear regression : y = 0.9126x + 15 mg/dL (0.15 g/L)

#### ❖ Limitations/Interférences (5-6)

- Do not report results outside of the usable range.

- According to SFBC recommendations some studies have been performed to determine the level of interference from different compounds :

**Bilirubin :** No significant interference up to 20 mg/dL (342,1 µmol/L).

**Triglycerides :** No significant interference up to 2500 mg/dL (28,25 mmol/L).

**Heparin :** No significant interference up to 50 mg/dL.

**Sodium citrate :** No significant interference up to 1000 mg/dL.

**Haemoglobin :** Negative bias from 0.8 g/L. As the physiological function of haptoglobin is to bind free hemoglobin, addition of hemoglobin leads to the formation of Hp-Hb complexes. This leads to falsely decreased results by 15% to 40%.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.<sup>(7)</sup>

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(8,9)</sup>

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

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**BIBLIOGRAPHIE/BIBLIOGRAPHY**

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1. Craig, W.Y. *et al.*, *Plasma Proteins: Clinical Utility and Interpretation*, Foundation for Blood Research, (2001), 126.
2. Johnson, A.M. *et al.*, *Proteins*, Tietz Fundamentals of Clinical Chemistry, 5<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 325.
3. Tietz, N.W., *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 306.
4. Dati, F., *et al.*, *Eur. J Clin Chem Clin Biochem*, (1996), 34, 517.
5. Vassault, A., *et al.*, *Ann. Biol. Clin.*, (1986), 44, 686.
6. Vassault, A., *et al.*, *Ann. Biol. Clin.*, (1999), 57, 685.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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