- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.<sup>(10)</sup>

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(11,12)</sup>

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

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### BIBLIOGRAPHIE/BIBLIOGRAPHY

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- Feldkamp, C.S., *Immunological reactions. Clinical Chemistry. Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A., Pe-scse, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 216.
- Craig, W.Y. *et al.*, *Plasma Proteins: Clinical Utility and Interpretation*. Foundation for Blood Research, (2001), 128.
- Johnson, A.M. *et al.*, *Proteins, Tietz Fundamentals of Clinical Chemistry*, 5<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 325.
- Tietz, N.W., *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 354.
- Dati, F., *et al.*, *Eur. J Clin Chem Clin Biochem*, (1996), 34, 517.
- Attalmanan, M., Levinson, S.S., *Clin. Chem.*, (2000), 46, 1230.
- Pontet, F., *et al.*, *Ann. Biol. Clin.*, (1997), 55, 486.
- Vassault, A., *et al.*, *Ann. Biol. Clin.*, (1986), 44, 686.
- Vassault A., *et al.*, *Ann. Biol. Clin.*, (1999), 57, 685.
- Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, *Acta Clin Belg.*, (2004), 59, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup>Ed., AACCPress, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup>Ed., AACCPress, (1995).

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