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BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Dufour, D.R., *The Pancreas: Function and Chemical Pathology, Clinical Chemistry, Theory, Analysis, Correlation*, 5th Ed., Kopian,L.A., Pesce, A.J., (Mosby, Inc.), (2010), 651, and appendix .
- Panteghini, M., Bais, R., *Enzyme, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
- Tietz NW et al. *Lipase in serum - the elusive enzyme: An overview. Clin Chem* (1993); **39**:746-756.
- Neumann U et al. *New substrates for the optical determination of lipase. EP 207252* (1987).
- Panteghini, M., *The never-ending search of an acceptable compromise for pancreatic lipase standardisation. Clin Chem Lab Med* 2012 ; **50**(3) : 419-421
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples.* (2002). WHO/DIL/ LAB/99.1 Rev.2.
- Wu, H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 634.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline.* CLSI (NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline.* CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCC Press, (1995).
- Glick MR, Ryder KW, Jackson SA. *Graphical comparisons of Interferences in Clinical Chemistry instrumentation.* *Clin Chem* 1986;**32**:471-475.

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English - EN

INTENDED USE

ELITech Clinical Systems LIPASE SL is an *in vitro* diagnostic reagent intended for the quantitative determination of lipase in human serum and plasma samples.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾


Lipase is a digestive enzyme of 48 kDa released by the pancreas which catalyses the hydrolysis of glycerol esters from triglycerides to form a monoglyceride and free fatty acids chains. Analysis of the activity of the lipase is mainly used in the diagnosis of the pancreatic disease (acute or chronic pancreatitis and their complication, carcinoma). During acute pancreatitis, a transitory increase of the activity of lipase is observed after 4 to 8h, reaches a peak after 24h, the activity returning normal after 8 to 14 days. However, an increase of the activity of lipase is also observed in other intra-abdominal pathologies: acute cholecystitis, pancreatic duct obstruction. Patients with a reduced glomerular filtration rate have also a increased lipase activity.

METHOD ⁽⁴⁾

Substrate: 1,2-O-Dilauryl-rac-Glycero-3-Glutaric acid (6-methylresorufin) ester (DGGM).
 Colorimetric- Kinetic.

PRINCIPLE ^(4,5)

The method for the determination of lipase is based on the cleavage of specific chromogenic lipase substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin) ester emulsified in stabilized micro-particles. In the presence of specific activators of pancreatic lipase as colipase, calcium ions and bile acids, the substrate is converted to 1,2-O-dilauryl-rac-glycerol and glutaric acid-(6-methylresorufin)ester which decomposes spontaneously to glutaric acid and methylresorufin. The increase of absorbance at 578 nm, due to methylresorufin formation, is proportional to the activity of lipase in the sample.

COMPOSITION

Reagent 1: R1

BICIN*buffer, pH 8.0	50 mmol/L
Colipase	≥ 0.9 mg/L
Sodium deoxycholate	1.6 mmol/L
Calcium chloride	10 mmol/L
Sodium azide	< 0.1 %

Détergent

Reagent 2: R2

Tartrate buffer, pH 4.16 (± 0.15)	10 mmol/L
1,2-O-dilauryl-rac-glycero-3-glutaric acid (6-methylresorufin)ester	0.27 mmol/L
Taurodeoxycholate	8.8 mmol/L

Détergent
 Preservative

*BICIN= N,N-bis(2-hydroxyethyl)glycine

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

WARNINGS AND PRECAUTIONS

- This *in vitro* diagnostic device is for professional use only.
- The reagent R1 contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.
- For more information, refer to the Safety Data Sheet (SDS).

STABILITIES

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES (6)

- Serum free from hemolysis
- Lithium heparinized plasma free from hemolysis.
- Do not use other specimens.

Warnings and precautions

According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.

Storage and stability

The samples are stable 7 days at room temperature, 3 weeks at 2-8 °C and 1 year at -20 °C.

REFERENCE VALUES (7)

Serum, plasma (37 °C): 13-60 U/L (0,22 - 1,00 µkat/L)

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

For ELITech Clinical Systems Selectra Analyzers

applications are available on request
 Wavelength 578 nm
 Temperature: 37 °C
 Read against reagent blank.

Reagent R1	150 µL
Sample	3 µL

Mix and wait 4 minutes 43 secondes incubation.

Reagent R2	90 µL
-------------------	-------

Mix and after 77 seconds incubation, measure the change of absorbance per minute (ΔA/min) during 159 seconds.

With Selectra TouchPro Software, use the application included in the barcode available at the end of this insert.

In the application, the offset must be set to: -2 U/L (-0.03 µkat/L)

Lipase SL reagent is strongly contaminated by Triglycerides SL reagent.

In order to avoid cuvette contamination on Selectra Pro instruments, program the following incompatibilities:

Software	Menu	Parameter
TouchPro	Test incompatibilities	Link / Triglycerides SL - Acid Solution
Other	Cuvette incompatibility	Triglycerides SL << HCl

In order to avoid needle contamination on Selectra Pro instruments, do not program Lipase SL and Triglycerides SL in the same run. Ensure the instrument goes back to "stand-by" status before launching a run containing Lipase SL.

CALCULATION

(ΔA) Sample x n n = calibrator concentration

(ΔA) Calibrator

Conversion factor : U/L x 0.0167 = µkat/L

CALIBRATION

For calibration, multiparametric calibrator ELICAL 2 must be used. Its value is traceable to the manual measurement.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

Measuring range

Determined according to CLSI® EP6-A protocol, the measuring range is from 5 to 300 U/L (0,08 to 5,00 µkat/L).

Samples having greater concentrations should be diluted 1:10 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 3000 U/L (50.00 µkat/L).

For Selectra TouchPro users, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

Determined according to CLSI EP17-A (9) protocol, LoD = 1 U/L (0,02 µkat/L).
 LoQ = 5 U/L (0,08 µkat/L).

Precision

Determined according to CLSI EP5-A2 protocol(10).

		Mean	Within-run	Total
	n	U/L	µkat/L	CV (%)
Level 1	80	28	0.47	0.8
Level 2	80	55	0.92	0.7
Level 3	80	229	3.82	0.4

Correlation

A comparative study has been performed between ELITEch Clinical systems method on ProM and another competitor system (Colorimetric method) on 100 human samples according to CLSI EP9-A2(11) protocol. The sample activities were between 6 and 284 U/L (0,10 to 4,73 µkat/L). The parameters of linear regression are as follows :
 Correlation coefficient : (r) = 0,999
 Linear regression : y = 1.063 x + 1 U/L (0,02 µkat/L)

Limitations/Interferences

- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according to EP7-A2(12) protocol of CLSI. Recovery within ± 10 % of initial value of lipase of 30, 60 and 240 U/L.

Unconjugated Bilirubin: No significant interference up to 30,0 mg/dL (513.1 µmol/L).
Conjugated Bilirubin: No significant interference up to 29.5 mg/dL (504.6 µmol/L).
Hemoglobin: No significant interference up to 50 mg/dL.

Triglycerides: No significant interference up to 3000 mg/dL (33.90 mmol/L).
Ascorbic acid: No significant interference up to 20 mg/dL.
Acetylsalicylic acid: No significant interference up to 200 mg/dL.
Acetaminophen: No significant interference up to 30 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.(13)

- Many other substances and drugs may interfere. Some of them are listed in Young and Glick. (14-16)

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

On board stability/Calibration frequency

On Board Stability: 14 days
Calibration frequency: 7 days
 Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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