PIT-MAL-3-10 (06/2018)

Français - FR

USAGE PRÉVU

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English - GB

INTENDED USE

ELITech Clinical Systems μALBUMIN IP is intended for the quantitative *in vitro* diagnostic determination of albumin in human urine samples.

CLINICAL SIGNIFICANCE (1-2)

Urinary excretion of albumin is generally significant of a renal damage. While in case of declared nephropathy the excretion of albumin is usually higher than 300 mg/day, an excretion rate ranging from 30 to 300 mg/day (called microalbuminuria) is often an early marker of nephropathy. Monitoring urinary albumin is thus recommended for diabetic individuals who are at high risk of developing nephropathies. Urinary albumin is also increased in case of physical exercise or during urinary infection.

METHOD

Immuno-turbidimetry.
End-Point.

PRINCIPE

The formation of albumin/anti-albumin antibody complexes, started by the addition of the antiserum to the sample, is accelerated in presence of polyethylene glycol. These complexes precipitate leading to an increase of turbidity measured at 340 nm. The albumin concentration is determined by means of a non linear calibration curve.

COMPOSITION

Reagent 1 : R1
Buffer pH 7,00
Polyethylene glycol 26 g/L
Reagent 2 : R2
Buffer pH 7,43
Polyclonal anti-human Albumin antibody (goat)

MATERIAL REQUIRED BUT NOT PROVIDED

- NaCl solution 9 g/L.
- Calibrators and Controls :

References	Description
IMAL-0043	μALBUMIN IP Calibrator Set
IMAL-0046	μALBUMIN IP Control I
IMAL-0047	μALBUMIN IP Control II

These references are sold separately

- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

PRECAUTIONS AND WARNING

- This reagent is for professional *in vitro* diagnostic use only.
- The reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

STABILITIES

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.
On board stability :
The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES (1,3)

Specimen
- Urine. Sample collection must not be done neither after muscular exercise nor in case of urinary tract infection. Samples which contain a precipitate must be centrifuged before analysis.
- A preliminary rapid test using urine strips is recommended to determine the concentration domain of the sample. If the concentration is higher than the analytical range, dilute the sample in a NaCl 9 g/L solution. The result must be multiplied by the dilution factor.
- Do not use other specimens.

Warnings and precautions

According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.
Storage and stability
Samples are stable for 7 days at 20-25 °C, 1 month at 2-8 °C or 6 months at -20 °C, with no added stabilizer.

REFERENCE VALUES (1)

Normal < 30 mg/day
Microalbuminuria 30-300 mg/day
Macroalbuminuria > 300 mg/day

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

This reagent can be used on most analysers, semi-automated analysers and manual method.
The applications are available on request.
Wavelength : 340 nm
Temperature : 37 °C
No dilution of samples

	CALIBRATION	TEST
Reagent R1	250 μL	250 μL
Calibrator	16 μL	-
Sample	-	16 μL

Mix and read the absorbances (A1) after 4 minutes and 40 seconds of incubation (sample blank) then add :

Reagent R2	40 μL
------------	-------

Mix and read the absorbance (A2) after 4 minutes incubation.

Note:

The results of μAlbumin urine samples may be falsely elevated due to contamination of the test by a serum sample. In order to avoid contamination on Selectra analysers, it is recommended to :

- calibrate the μAlbumin IP reagent separately from calibration of seric tests.
- program the μAlbumin IP test in a run that does not contain any serum test. Alternatively re-measure only the pathological results in a separate run.

Ensure the instrument goes back to "stand-by" status before launching any separate run containing μAlbumin IP tests.

CALCULATION

Albumin concentration is calculated from a calibration curve obtained from five calibrators of different levels and a zero point.
Conc = f (ΔA) = f (A2 - A1)

CALIBRATION

For calibration, use ELITech μAlbumin IP Calibrator set, ready-to-use.
Add a zero point (NaCl 9 g/L solution)

Concentration values for calibrators are traceable to ERM-DA470k/IFCC reference material .

It is recommended to make a new calibration when reagent lots change and when quality control results fall outside the established range.

QUALITY CONTROL

To check the accuracy of assays, control sera such as μAlbumin IP Control I and Control II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.

PERFORMANCE DATA at 37 °C on Selectra Junior

- Analytical range
The reagent is linear from 0.30 to 40,00 mg/dL (3.0 to 400.0 mg/L).
The exact range depends on the value of the calibrator used.

- Detection limit (5)
Determined according to SFBC protocol, the detection limit is equal to 0.08 mg/dL (0.8 mg/L).

...

- Hook effect
No risk of error linked to hook effect up to 200 mg/dL (2000 mg/L).

- Precision

Within-run reproducibility	n	Mean		CV (%)
		mg/dL	mg/L	
Human urine N° 1	20	2.50	25.0	2.6
Human urine N° 2	20	5.04	50.4	1.3
Human urine N° 3	20	21.60	216.0	2.3

Between-run reproducibility	n	Mean		CV (%)
		mg/dL	mg/L	
Human urine N° 1	20	2.76	27.6	4.1
Human urine N° 2	20	5.54	55.4	4.9
Human urine N° 3	20	21.94	219.4	2.6

- Correlation
A comparative study was performed on 33 human urine samples. The values were between 0,49 and 34,29 mg/dL (4,9 and 342,9 mg/L).
The parameters of linear regression are as follows :
Correlation coefficient : (r) = 0,9817
Linear regression : y = 0,9119 x - 0,52 mg/dL (5,2 mg/L)

- Limitations/Interferences (5)
Do not report results outside of the usable range.

- According to SFBC recommendations some studies have been performed to determine the level of interference from different compounds :

Conjugated Bilirubin : No significant interference up to 23,4 mg/dL (400 μmol/L).
Unconjugated Bilirubin : No significant interference up to 36 mg/dL (616 μmol/L).
Haemoglobin : No significant interference up to 500 mg/dL (5 g/L).

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. (6-7)

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- On board stability/Calibration frequency (4-5)
On Board Stability: 14 days
Calibration frequency: 14 days

Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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BIBLIOGRAPHIE/BIBLIOGRAPHY/ BIBLIOGRAFÍA

1. Sacks, D.B., *Carbohydrates, Tietz Fundamentals of Clinical Chemistry*, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 427.
2. First, M.R., *Renal function. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 477 and appendix.
3. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*.(2002). WHO/DIL/LAB/99.1 Rev.2.
4. Vassault, A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
5. Vassault A., et al., *Ann. Biol. Clin.*, (1999), **57**, 685.
6. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2ndEd., AACC Press, (1997).
7. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4thEd., AACC Press, (1995).

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