VTL-PRTU-4-v9 (07/2019)

**Français - FR**

☞ **USAGE PRÉVU**

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**English - EN**

☞ **INTENDED USE**

ELITech Clinical Systems MICROPROTEIN PLUS is an *in vitro* diagnostic reagent intended for the quantitative determination of total protein in human urine samples, using MICROPROTEIN PLUS Standard 100 mg/dL.

**CLINICAL SIGNIFICANCE** <sup>(1-3)</sup>

An increase of urinary protein (normally very low) most often indicates a renal disorder due to an increased glomerular permeability (nephrotic syndrome, progressive renal failure, hypertension) or to a decreased tubular reabsorption (Fanconi syndrome, chronic pyelonephritis, heavy metal poisoning). There are also pathological proteinurias due to protein overload or to disorders of the lower urinary tract. Furthermore, an increased excretion of albumin for diabetics of type 1 is an early marker of renal disease. Under some conditions (intense exercise, strong fever, prolonged upright posture, pregnancy...), the urinary protein rate can be increased in an intermittent way.

**METHOD** <sup>(4)</sup>  
Colorimetric: Pyrogallol Red.  
End point.

**PRINCIPLE** <sup>(4)</sup>

Proteins bound to the pyrogallol red-molybdate complex form in an acidic medium a blue/purple complex whose absorbance is proportional to the protein concentration in the sample.

☞ **COMPOSITION**

**Reagent : R**

Succinic acid, pH 2.5	50	mmol/L
Pyrogallol red	60	µmol/L
Sodium molybdate	165	µmol/L
Sodium oxalate	1	mmol/L

**Standard : Std** (ref: PRTU-0600)

Albumin	100	mg/dL
	1000	mg/L
Sodium azide	< 0.1	%

**MATERIALS REQUIRED BUT NOT PROVIDED**

- PRTU-0022, MICROPROTEIN PLUS Standard 100 mg/dL for the kit Ref PRTU-0250.
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

☞ **WARNINGS AND PRECAUTIONS**

- These *in vitro* diagnostic devices (reagent and standard) are for professional use only.  
- The reagent R is classified as hazardous (methanol):

**WARNING** : May cause damage to organs. Do not breathe mist/vapours/spray. Wash thoroughly after handling. IF exposed or concerned: Call a POISON CENTER/doctor.

- The standard contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of this standard always flush with copious amounts of water to prevent azide buildup.
- For more information, refer to the Safety Data Sheet (SDS).
- Take normal precautions and adhere to good laboratory practice.
- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- Use clean or single use laboratory equipment only to avoid contamination.

☞ **STABILITIES**

**Store at 2-8 °C and protect from light. Do not freeze.**  
Do not use after expiration dates indicated on the vial labels.  
**On board stability** : The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

**PREPARATION**

The reagent and standard are ready to use.

☞ **PRODUCT DETERIORATION**

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

**SAMPLES** <sup>(5)</sup>

**Specimen**  
- Urine, 24h collection recommended.  
- Do not use other specimens.

**Warnings and precautions**  
- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.

- Do not use samples containing blood.
- Storage**  
- During 24h collection, urine should be stored at 2-8 °C. Then, urine should be analysed fresh.  
When storage is needed, urine should be stored :
  - 1 day at room temperature.
  - 7 days at 2-8 °C.
  - 1 month at -20 °C.

**REFERENCE VALUES** <sup>(1,3)</sup>

Urine : 1-14 mg/dL (10-140 mg/L).  
The values are strongly increased after an intense muscular exercise.

*Note* : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.