☛ **Références/Referencen/ Referencias/ Referências:** **Compositi** **on du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :**
 PROB-0250 R 12 x 20 mL
 PROB-0600 R 2 x 125 mL + Std 1 x 5 mL
 PROB-0700 R 4 x 250 mL + Std 1 x 5 mL

☛ **Références/Referencen/ Referencias/ Referências:** **Compositi** **on du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :**
 PROB-0250 R 12 x 20 mL
 PROB-0600 R 2 x 125 mL + Std 1 x 5 mL
 PROB-0700 R 4 x 250 mL + Std 1 x 5 mL

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**BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFIA/BIBLIOGRAFIA**

- Scherwin, J.E., *Liver function. Clinical Chemistry. Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
- Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 518.
- Doumas, B.T., *et al., Clin. Chem.*, (1981), **27**, 1642.
- Guder, W.G., *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Vassault, A., *et al., Ann. Biol. Clin.*, (1986), **44**, 686.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

☛ **SYMBOLES/SYMBOLS/ SÍMBOLOS/SIMBOLOS**
 - Les symboles utilisés sont décrits dans la norme ISO-15223-1 hormis ceux présentés ci-dessous.
 - Symbols used are defined on ISO-15223-1 standard, except those presented below.
 - Los símbolos utilizados son descritos en la norma ISO-15223-1 a la excepción de los presentados a continuación.
 - Os símbolos utilizados são definidos na norma ISO-15223-1, exceto os apresentados abaixo.

CONT	Contient Content Contiene Conteúdo
R	Réactif Reagent Reactivo Reagente
Std	Standard Standard Estándar Padrão
CE	Conformité Européenne European Conformity Conformidad Europea Conformidade Europeia

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IMPORTANT NOTE
 - Only for ref. PROB-0250, used with Selectra TouchPro software.
 - See § PROCEDURE: Contamination risk

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English - EN

☛ **INTENDED USE**

ELITech Clinical Systems TOTAL PROTEIN PLUS is an *in vitro* diagnostic reagent intended for the quantitative determination of total protein in human serum and plasma samples.

CLINICAL SIGNIFICANCE ⁽¹⁻²⁾

In human plasma, albumin accounts for 50 to 60 % of total proteins ; the reminder fraction mainly contains globulins (α 1, α 2, β and γ). Most plasmatic proteins are synthesized by the liver, except immunoglobulins.

Increase of the plasmatic volume (salt retention syndrome, intoxication with water...) or its reduction (dehydration related to vomiting, diarrhoea...) induce respectively relative hypoproteinemina and relative hyperproteinemia.

For a normal plasmatic volume, abnormal total protein rates only occur in the event of disorder affecting the concentration of albumin or immunoglobulins. Thus, severe proteinic insufficiency (malabsorption, maldigestion, dietary insufficiency), renal and hepatic diseases result in hypoproteinemia. If total protein concentration is lower than 4 g/dL, oedemas can be observed. Hyperproteinemia can be seen, for example, in case of hyperimmunoglobulinemia (multiple myeloma, infection).

METHOD ⁽³⁾

Biuret. End point

PRINCIPLE ⁽³⁾

Serum proteins form a coloured complex in the presence of copper salt in alkaline solution.

- For more information, refer to the Safety Data Sheet (SDS).
- The standard contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of this standard always flush with copious amounts of water to prevent azide buildup.
- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.

➤STABILITIES
Store at 2-25 °C and protect from light. Do not freeze.

Do not use after expiration dates indicated on the vial labels.
On board stability :
The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION
The reagent and standard are ready to use.**➤PRODUCT DETERIORATION**
- The Product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).**SAMPLES** ^(1,2,4)
Specimen
- Serum.
- Lithium heparinized plasma.
- Do not use other specimens.
➤Warnings and precautions
- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.
- Samples must be free from haemolysis and lipemia.
Storage and stability
Samples are stable for 7 days at 2-8 °C and at least 2 months at -20 °C. For longer storage, freeze samples at -70 °C.**REFERENCE VALUES** ^(1,2,4)
Serum: Ambulatory patients 6,4 - 8,3 g/dL
Patients at rest 6,0 - 7,8 g/dL
64 - 83 g/L 60 - 78 g/L

Plasma:
Plasma concentrations are increased from 0.2 to 0.4 g/dL (2 to 4 g/L) compared to serum concentrations (fibrinogen).

Note : *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

➤PROCEDURE
Manual Procedure
Wavelength : 546 nm
Optical path : 1 cm
Sample/reagent ratio : 1:100
Temperature : 37 °C
Read against reagent blank.

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/ Calibrator	10 µL	-
Sample	-	10 µL

Mix and read the absorbances (A) after an incubation of 10 minutes.

Automatic Procedure
These reagents may be used in several automatic analyzers. For ELITech Selectra ProM Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.**Important set-up information:**
MAGNESIUM XYLIDYL reagent can be weakly contaminated by TOTAL PROTEIN PLUS on Selectra ProM and ProXL.

In order to avoid contamination on these instruments , program the following incompatibilities:

Software	Menu	Parameter
TouchPro	Probe incompatibilities	incompatibility/ PROTEIN - MAGNESIUM
Other	Needle incompatibility	PROTEIN : MAGNESIUM

CALCULATION
$$\frac{A_{\text{Sample}}}{A_{\text{Standard/ Calibrator}}} \times n \quad n = \text{calibrator/standard concentration}$$

Conversion factor : g/dL x 10 = g/L**CALIBRATION**
For the references PROB-0600/PROB-0700 : For calibration, use either multiparametric calibrator ELICAL 2 or Total Protein Standard 6 g/dL.
For the reference PROB-0250 : For calibration , use multiparametric calibrator ELICAL 2.

Concentration values of Total Protein Standard 6 g/dL and multiparametric calibrator ELICAL 2 are traceable to the Standard Reference Material SRM 909c (of the National Institute of Standards and Technology).

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).**➤QUALITY CONTROL**
To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.**WASTE MANAGEMENT**
Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.**PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers**
- Measuring range
Determined according to CLSI EP6-A protocol⁽⁵⁾ the measuring range is from 0.20 to 12.0 g/dL (2.0 to 120.0 g/L).**- Limit of Detection (LoD) and Limit of Quantification (LoQ)**
Determined according to CLSI EP17-A protocol⁽⁶⁾.
LoD = 0.03 g/dL (0.3 g/L)
LoQ = 0.10 g/dL (1.0 g/L)**- Precision**
Determined according to CLSI EP5-A2 protocol⁽⁷⁾.

		Mean	Within-run	Total	
	n	g/dL	g/L	CV (%)	
Low level	80	4.03	40.3	0.4	1.0
Medium level	80	6.62	66.2	0.3	1.6
High level	80	9.06	90.6	0.5	1.1

- Correlation
A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (Biuret method) on 100 human sera samples according to CLSI EP9-A2 protocol⁽⁸⁾. The values were between 0.27 and 11.25 g/dL (2.7 and 112.5 g/L). The parameters of the linear regressions are as follows:
Correlation coefficient : (r) = 0.997
Linear regression: y = 0.993 x + 0.05 g/dL (0.5 g/L)**- Limitations/Interferences**
- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol⁽⁹⁾ and SFBC recommendations⁽¹⁰⁾. Recovery is within ±10% of initial value at total protein concentration of 4.00; 6.50 and 9.00 g/dL.

Unconjugated Bilirubin : No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated Bilirubin: No significant interference up to 29.5 mg/dL (504 µmol/L).
Glucose: No significant interference up to 507 mg/dL (28.14 mmol/L).
Turbidity: No significant interference up to 263 mg/dL (2.97 mmol/L) triglycerides equivalent.
Hemoglobin: No significant interference up to 300 mg/dL.
Dextran: Induces falsely high results at therapeutic concentrations.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽¹¹⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽¹²⁻¹³⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- On board stability/Calibration frequency
On Board Stability: 14 days
Calibration frequency: 14 days
Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.**➤These performances have been obtained using ELITech Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.**

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