

Consult Instruction for Use.

APPLICATION COBAS MIRA - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

GENERAL

Measurement Mode : ABSORB
 Reaction Mode : R-S
 Calibration Mode : CALIBRATOR
 Reagent Blank : REAG/DIL
 Cleaner : NO
 Wavelength : 500 nm
 Decimal position : 0
 Unit : mg/dL

ANALYSIS

Sample Dil.Name : H₂O
 Post Dil.Factor : NO
 Conc.Factor : NO
Sample Cycle 1
 Vol : 3 µL Dil : 10 µL
Reagent Cycle 1
 Vol : 290 µL

CALCULATION

Sample limit : NO
 Reac.Direction : INCREASE
 Check : ON
 Convers.Factor : 1.00000
 Offset : 0.00000
 Test Range Low : *
 High : *
 Norm.Range Low : *
 High : *
 Number of steps : 1
 Calc.Step A : KINETIC
 Readings First : 1
 Last : 4
 Reaction Limit : NO

CALIBRATION

Calib.Interval : ON REQUEST
 BLANK
 Reag.Range Low : -0.150***
 High : 0.450***
 Blank Range Low : -0.100
 High : 0.150
 STANDARD POS : As desired
 VAL : **
 REPL : TRIPL
 DEV : 5%

* : Entered by the operator.
 ** : Entered standard value
 ***: Defined for a concentration of 100 mg/dL
 These data must be modified for other concentration.