

#### **ELITech Clinical Systems SAS**

# **CHOLESTEROL LDL SL 2G**

REF: LDLL-XXXX

### <u>Instruction for use</u>:

Consult Instruction for Use

## APPLICATION - MINDRAY BS-120 / BS-200 / BS-200E - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

### BS120 / BS200 / BS200E

PARAMETERS			
Test	LDLL	R1	300 / 300 / 240
No.	**	R2	100 / 100 / 80
Full Name	LDL	Sample volume	3/3/2.4
Standard No.	1	R1 Blank	
Reac. Type	Endpoint	Mixed Rgt Blank	
Pri. Wavelength	578 nm / 578 nm / 570 nm	Linearity Range	*
Sec. Wavelength		Linearity Limit	*
Direction	Increase	Substrate Limit	*
Reac Time	0 20	Factor	
Incuba. Time	19	Prozone check	*
Units	mg/dL	q1	Q2
Precision	Integer	q3	Q4
	·	PC	Abs
CALIBRATION (Cal + Rgt Blk)			
Rule	Linear one point / Linear two point / Linear two point		
Sensitivity	1		
Replicates	2		
Interval (days) Difference I imit	0		
SD			
Blank Response			
Error Limit			
Correlation coefficient	1		

<sup>\*</sup>Enter by the user.