

# ALP (DEA) SL

## ONE REAGENT PROCEDURE

REF : PASL-XXXX

Instruction for use :  
 Consult Instruction for Use

### APPLICATION - MINDRAY BS-120 / BS-200 / BS-200E - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

**Working reagent preparation:** Mix 4 volumes of reagent R1 with 1 volume of reagent R2.  
**Working Reagent Stability:** 5 days at 20-25°C / 4 weeks at 2-8°C

### BS120 / BS200 / BS200E

PARAMETERS			
Test No.	PASL **	R1 R2	200 / 200 / 200
Full Name	ALP	Sample volume	4 / 4 / 4
Standard No.	1	R1 Blank	
Reac. Type	Kinetic	Mixed Rgt Blank	
Pri. Wavelength	405 nm / 405 nm / 412 nm	Linearity Range	*
Sec. Wavelength		Linearity Limit	*
Direction	Increase	Substrate Limit	*
Reac Time	3 8	Factor	
Incuba. Time		Prozone check	*
Units	U/L	q1	Q2
Precision	Integer	q3	Q4
		PC	Abs
CALIBRATION (Cal + Rgt Blk)			
Rule	Linear one point / Linear two point / Linear two point		
Sensitivity	1		
Replicates	2		
Interval (days)	0		
Difference Limit			
SD			
Blank Response			
Error Limit			
Correlation coefficient	1		

\*Enter by the user.