

Consult Instruction for Use

REF : CKSL-XXXX

APPLICATION - Mindray - BS 360 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Define / Edite Chemestries

Chem: CKSL N°: Sample Type: Serum

Chemistry: CK Print name:

Reaction Type: Kinetic Reaction Direction: Positive

Pri Wave: 340 nm Sec Wave:

Unit: U/L Decimal: 0

Incubation Time: 28

Blank time: Reaction Time: 13 11

Sample Vol Aspirated Diluent Reagent Vol

Standard 6 μL μL μL R1 120 μL

Decreased μL μL μL R2 30 μL

Increased μL μL μL

Sample Blank Auto rerun

Define / Edite Chemestries

Chem: CKSL N°: Sample Type: Serum

Chemistry: CK Print name:

Linearity range (standard) * Linearity Limit *

Linearity range (Decreased) * Substrate Depletion

Linearity range (Increased) * Mixed Blank Abs *

R1 Blank Abs * On-board Stability * Jour

Blank Response * Reagent Alarm Limit

Twin Chemistry Enzyme Linear Extension

Prozone Check Rtae Check

Q1 Q2 Q3 Q4 PC ABS

Calibration setup

Chem: CKSL

Calibration settings

Math Model: Two-point linear

Factor: Replicates: 2

Acceptance Limits

Cal Time: Heure

Slope Diff: SD:

Sensitivity: Reproductb.:

Deter Coeff:

Auto calib.

Bottle Changed Lot Changed Cal Time

Calibrator	Conc.	Pos	Lot No.
Water	0	W	
Cal1	xxx	xxx	

*Values entered by the operator