

Consult Instruction for Use

APPLICATION - Mindray - BS 360 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Define / Edite Chemestries

Chem	CLDL	N°		Sample Type	Serum
Chemistry	LDL			Print name	
Reaction Type	Endpoint	Reaction Direction	Positive		
Pri Wave	570 nm	Sec Wave	700 nm		
Unit	mg/dL	Decimal	0		
Blank time	-3	-1	Incubation Time	28	
			Reaction Time	39 41	
Sample Vol	Aspirated	Diluent	Reagent Vol		
Standard	2 μL		R1	150 μL	
Decreased			R2	50 μL	
Increased					
	<input type="checkbox"/> Sample Blank	<input type="checkbox"/> Auto rerun			

Define / Edite Chemestries

Chem	CLDL	N°		Sample Type	Serum
Chemistry	LDL			Print name	
Linearity range (standard)	*	*	Linearity Limit	*	
Linearity range (Decreased)	*	*	Substrate Depletion		
Linearity range (Increased)	*	*	Mixed Blank Abs	*	
R1 Blank Abs	*	*	On-board Stability	*	Jour
Blank Response	*	*	Reagent Alarm Limit		
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rtae Check				
Q1	Q2	Q3	Q4	PC	ABS

Calibration setup

Chem

Calibrator	Conc.	Pos	Lot No.
Water	0	W	
Cal1	xxx	xxx	

Calibration settings	Factor	Replicates
Math Model	<input type="text" value="Two-point linear"/>	<input type="text" value="2"/>
Acceptance Limits	Cal Time	Heure
Slope Diff	<input type="text"/>	SD <input type="text"/>
Sensitivity	<input type="text"/>	Reproductb. <input type="text"/>
Deter Coeff	<input type="text"/>	
Auto calib.	<input type="checkbox"/> Bottle Changed	<input type="checkbox"/> Lot Changed
	<input type="checkbox"/> Cal Time	

*Values entered by the operator