

APPLICATION - Mindray - BS 360 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Define / Edite Chemestries

Chem	PASL	N°		Sample Type	Serum
Chemistry	ALP			Print name	
Reaction Type	Kinetic	Reaction Direction	Positive		
Pri Wave	412 nm	Sec Wave			
Unit	U/L	Decimal	0		
		Incubation Time	28		
Blank time		Reaction Time	5	13	
Sample Vol		Aspirated		Diluent	
Standard	3 μ L			Reagent Vol	
Decreased				R1	120 μ L
Increased				R2	30 μ L
<input type="checkbox"/>	Sample Blank	<input type="checkbox"/>	Auto rerun		

Define / Edite Chemestries

Chem	PASL	N°		Sample Type	Serum
Chemistry	ALP			Print name	
Linearity range (standard)	* <input type="text"/>	* <input type="text"/>	Linearity Limit	* <input type="text"/>	
Linearity range (Decreased)	* <input type="text"/>	* <input type="text"/>	Substrate Depletion	<input type="text"/>	
Linearity range (Increased)	* <input type="text"/>	* <input type="text"/>	Mixed Blank Abs	* <input type="text"/>	<input type="text"/>
R1 Blank Abs	* <input type="text"/>	* <input type="text"/>	On-board Stability	* <input type="text"/>	Jour
Blank Response	* <input type="text"/>	* <input type="text"/>	Reagent Alarm Limit	<input type="text"/>	
Twin Chemistry	<input type="text"/>			<input type="checkbox"/>	Enzyme Linear Extension
<input type="checkbox"/>	Prozone Check	<input checked="" type="radio"/>	Rtae Check		
Q1	<input type="text"/>	Q2	<input type="text"/>	Q3	<input type="text"/>
Q4	<input type="text"/>	PC	<input type="text"/>	ABS	<input type="text"/>

Calibration setup

Chem	PASL
Calibration settings	
Math Model	Two-point linear
Factor	<input type="text"/>
Replicates	2
Acceptance Limits	
Cal Time	<input type="text"/> Heure
Slope Diff	<input type="text"/>
SD	<input type="text"/>
Sensitivity	<input type="text"/>
Reproductb.	<input type="text"/>
Deter Coeff	<input type="text"/>
Auto calib.	
<input type="checkbox"/>	Bottle Changed
<input type="checkbox"/>	Lot Changed
<input type="checkbox"/>	Cal Time

Calibrator	Conc.	Pos	Lot No.
Water	0	W	
Cal1	xxx	xxx	

*Values entered by the operator