

Instruction for use :  
 Consult Instruction for Use

### APPLICATION - PRESTIGE24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Item Name FERR <u>DATA INFORMATION</u> Units ug/L Decimals 0 <u>ANALYSIS</u> Type END W.Length 1 546 Method <u>CORR</u> SLOPE INTER 1.000 x + 0		<u>CALIBRATION</u> TYPE Spline STANDARD #1 0,0625 x Cal. Val #4 0,50 x Cal. Val #2 0,125 x Cal. Val #5 1,00 x Cal. Val #3 0,250 x Cal. Val #6 <u>NORMAL RANGE</u> LOW HIGH SERUM MALE FEMALE URINE	
Item Name FERR <u>ASPIRATION</u> KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 18 µL REAGENT 1 160 µL REAGENT 2 40 µL Third Mix <input checked="" type="checkbox"/> OF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B		<u>DATA PROCESS</u> <u>READ</u> START END MAIN 51 52 SUB 33 34 ENDPOINT LIMIT 3 LINEAR CHECK (%) <u>FACTOR</u> Blank Correction 1.000 <u>PROZONE CHECK</u> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	
<u>MONITOR</u> 0 LEVEL POINT 1 SPAN 3.000			

**\*\* Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.**