


ALBUMIN

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Instruction for use : 

Working temperature : 37°C

| | | | |
|-----------|--|-----|------|
| Item Name | | XXX | ALBU |
|-----------|--|-----|------|

| | |
|------------------|------|
| DATA INFORMATION | |
| UNITS | g/dL |
| DECIMALS | 1 |

| | |
|-----------------|-----|
| ANALYSIS | |
| TYPE | END |
| Main W Length 1 | 600 |
| Sub W Length 1 | |
| METHOD | |

| | | |
|------|-------|-------|
| CORR | SLOPE | INTER |
| Y= | 1 | X+ 0 |

| | |
|-------------|--------|
| CALIBRATION | |
| TYPE | Linear |
| #1 | xxx #4 |
| #2 | #5 |
| #3 | #6 |

| | | | | |
|--------------|------|------|--------|------|
| NORMAL RANGE | | | | |
| | MALE | | FEMALE | |
| | Low | High | Low | High |
| Serum | xxx | xxx | xxx | xxx |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| other | | | | |

| | | | |
|-----------|--|-----|------|
| Item Name | | XXX | ALBU |
|-----------|--|-----|------|

| | |
|---------------|---|
| Aspiration | |
| KIND | <input type="radio"/> single <input type="radio"/> double |
| SAMPLE | VOLUME |
| REAGENT 1 vol | 350 µL |

| | |
|--------------|-----------|
| DATA PROCESS | |
| READ | Start End |
| Main | 17 19 |
| Sub | |

| | |
|------------------|--------|
| ABSORBANCE LIMIT | |
| Low | -3.000 |
| High | 2.000 |

| | |
|------------------|---|
| FACTOR | |
| Blank correction | 1 |

| | | |
|------------------|--|---|
| Endpoint limit | | 2 |
| Linear check (%) | | 0 |

| | |
|-----------|--|
| Third mix | <input type="radio"/> OFF <input type="radio"/> ON |
|-----------|--|

| | |
|----------|--|
| R1 Blank | <input type="radio"/> Water-blank <input type="radio"/> R1-blank-1 |
|----------|--|

| | |
|----------|--|
| Dilution | |
| Diluent | <input type="radio"/> 99:Di1 <input type="radio"/> 100:Di2 |

| | | | |
|---------------|-------|-----|---|
| PROZONE CHECK | | | |
| | start | End | Limit % |
| First | | | low high |
| Second | | | <input type="radio"/> <input type="radio"/> |
| third | | | <input type="radio"/> <input type="radio"/> |

| | |
|---------------|---|
| Monitor | |
| 0 level point | 1 |
| SPAN | 3 |


| | | | |
|-----------|--|-----|------|
| Item name | | XXX | ALBU |
|-----------|--|-----|------|

| | | | |
|--|--|-----------------------------------|--|
| Auto Rerun SW | | Auto rerun condition (Absorbance) | |
| <input type="radio"/> ON <input type="radio"/> OFF | | ABS. range | |
| | | lower | <input type="radio"/> ON <input type="radio"/> OFF |
| | | high | <input type="radio"/> ON <input type="radio"/> OFF |

| | | | |
|--|--|----------------|--|
| Auto Rerun Range (result) | | Prozones range | |
| <input type="radio"/> ON <input type="radio"/> OFF | <input type="radio"/> ON <input type="radio"/> OFF | | <input type="radio"/> ON <input type="radio"/> OFF |

| | | |
|----------|-------|--------|
| | lower | higher |
| Serum | 1.5 | 6.0 |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| other | | |

XXX: enter data by the user

 Modification from the previous version