


ALT / GPT 4+1 SL

TWO REAGENT PROCEDURE

Ref. : ALSL-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Instruction for use : 

Working temperature : 37°C

Item Name XXX ALSL		CALIBRATION	
DATA INFORMATION		TYPE Factor	
UNITS	L/L	Factor	-1746
DECIMALS	0		
ANALYSIS		NORMAL RANGE	
TYPE RATE		MALE FEMALE	
Main W Length 1	340	Low High	Low High
Sub W Length 2		Serum XXX XXX	XXX XXX
METHOD		Urine	
		Plasma	
		CSF	
		Dialysis	
		other	
CORR	SLOPE	INTER	
Y=	1	X+	0

Item Name XXX ALSL		DATA PROCESS		ASPIRATION	
KIND <input type="radio"/> single <input checked="" type="radio"/> double		READ		LIMIT	
SAMPLE	VOLUME	Main	Start End	Low	0.800
REAGENT 1 vol	200	Sub	36 45	High	3.000
REAGENT 2 vol	50	FACTOR		Endpoint limit	0
		Blank correction 1		Linear check (%)	40
Third mix	<input type="radio"/> OFF <input type="radio"/> ON	Dilution			
R1 Blank	<input type="radio"/> Water-blank <input type="radio"/> R1-blank-1	Diluent <input type="radio"/> 99-Dil 1 <input type="radio"/> 100-Dil 2			
Monitor		PROZONE CHECK			
0 level point	1	start	End	Limit %	
SPAN	3	First			low high
		Second			<input type="radio"/> <input type="radio"/>
		third			<input type="radio"/> <input type="radio"/>

Item name XXX ALSL		Auto rerun condition (Absorbance)	
Auto Rerun SW <input checked="" type="radio"/> ON <input type="radio"/> OFF		ABS. range	
Auto Rerun Range (result) <input checked="" type="radio"/> ON <input type="radio"/> OFF		lower	<input type="radio"/> ON <input type="radio"/> OFF
		high	<input type="radio"/> ON <input type="radio"/> OFF
		Prozone range	<input type="radio"/> ON <input type="radio"/> OFF

XXX: enter data by the user