


AMYLASE SL

Réf. : AMSL-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	AMSL
-----------	--	-----	------

DATA INFORMATION

UNITS	U/L
DECIMALS	0

ANALYSIS

TYPE	RATE
Main W Length 1	405
Sub W Length 1	

METHOD

METHOD	
--------	--

CORR

SLOPE	INTER
Y= 1	X+ 0

CALIBRATION

TYPE	Factor
Factor	4640

NORMAL RANGE

	MALE		FEMALE	
	Low	High	Low	High
Serum	XXX	XXX	XXX	XXX
Urine				
Plasma				
CSF				
Dialysis				
other				

Item Name		XXX	AMSL
-----------	--	-----	------

Aspiration

KIND single double

SAMPLE	VOLUME	
REAGENT 1 vol	200	µL

DATA PROCESS

READ	Start	End
Main	17	29
Sub		

ABSORBANCE LIMIT

Low	-3.000
High	2.000

FACTOR

Blank correction	1
------------------	---

Endpoint limit

Endpoint limit	0
Linear check (%)	40

Third mix OFF ON

R1 Blank Water-blank R1-blank-1

Dilution

Diluent 99-Dil 1 100-Dil 2

PROZONE CHECK

	start	End	Limit %	low	high
First				<input type="radio"/>	<input type="radio"/>
Second				<input type="radio"/>	<input type="radio"/>
third				<input type="radio"/>	<input type="radio"/>

Monitor

0 level point	1
SPAN	3

Item name		XXX	AMSL
-----------	--	-----	------

Auto Rerun SW ON OFF

Auto Rerun Range (result) ON OFF

Auto rerun condition (Absorbance)

ABS. range lower ON OFF

high ON OFF

Prozone range ON OFF

	lower	higher
Serum	25	1600
Urine		
Plasma		
CSF		
Dialysis		
other		

XXX: enter data by the user