


AST / GOT 4+1 SL

TWO REAGENT PROCEDURE

Ref. : ASSL-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	ASSL
-----------	--	-----	------

DATA INFORMATION
 UNITS: U/L
 DECIMALS: 0

ANALYSIS
 TYPE: RATE
 Main W.Length 1: 340
 Sub W.Length 2:
 METHOD:

CORR: SLOPE: 1 INTER: X+ 0

CALIBRATION
 TYPE: Factor
 Factor: -1746

	MALE		FEMALE	
	Low	High	Low	High
Serum	XXX	XXX	XXX	XXX
Urine				
Plasma				
CSF				
Dialysis				
other				

Item Name		XXX	ASSL
-----------	--	-----	------

Aspiration
 KIND: single double

DATA PROCESS
 READ: Start End
 Main: 36 45
 Sub:
 REAGENT 1 vol: 200 µL
 REAGENT 2 vol: 50 µL

ABSORBANCE LIMIT
 Low: 0.600
 High: 3.000

FACTOR
 Blank correction: 1

Endpoint limit: 0
 Linear check (%): 40

Third mix: OFF ON

R1 Blank: Water-blank R1-blank-1

Dilution
 Diluent: 99: Dil 1 100: Dil 2

PROZONE CHECK

	start	End	Limit %	low	high
First				<input type="radio"/>	<input type="radio"/>
Second				<input type="radio"/>	<input type="radio"/>
third				<input type="radio"/>	<input type="radio"/>

Monitor
 0 level point: 1
 SPAN: 3

Item name		XXX	ASSL
-----------	--	-----	------

Auto Rerun SW: ON OFF

Auto Rerun Range (result): ON OFF

Auto rerun condition (Absorbance)
 ABS. range: lower high
 Prozone range: ON OFF

	lower	higher
Serum	10	250
Urine		
Plasma		
CSF		
Dialysis		
other		

XXX: enter data by the user