


# URIC ACID MONO SL

Réf. : AUML-XXXX

## APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

**Instruction for use** : 

**Working temperature** : 37°C

Item Name		XXX	AUML
-----------	--	-----	------

**DATA INFORMATION**

UNITS	ma/dL
DECIMALS	2

**ANALYSIS**

TYPE	END
Main W Length 1	546
Sub W Length 1	

**METHOD**

METHOD	
--------	--

**CORR**

SLOPE	INTER
Y= 1	X+ 0

**CALIBRATION**

TYPE	Linear		
#1	xxx	#4	
#2		#5	
#3		#6	

**NORMAL RANGE**

	MALE		FEMALE	
	Low	High	Low	High
Serum	xxx	xxx	xxx	xxx
Urine				
Plasma				
CSF				
Dialysis				
other				

---

Item Name		XXX	AUML
-----------	--	-----	------

**ASPIRATION**

KIND	<input type="radio"/> single <input type="radio"/> double
SAMPLE	VOLUME 5 µL
REAGENT 1 vol	200

**DATA PROCESS**

READ	Start	End
Main	35	38
Sub		

**ABSORBANCE LIMIT**

Low	-3.000
High	2.00

**FACTOR**

Blank correction	1
Endpoint limit	2
Linear check (%)	0

**Third mix**  OFF  ON

**R1 Blank**  Water-blank  R1-blank-1

**Dilution**

Diluent	<input type="radio"/> 99:Di1 <input type="radio"/> 100:Di2
---------	--

**PROZONE CHECK**

	start	End	Limit %	low	high
First				<input type="radio"/>	<input type="radio"/>
Second				<input type="radio"/>	<input type="radio"/>
third				<input type="radio"/>	<input type="radio"/>

**Monitor**

0 level point	1
SPAN	3

---

Item name		XXX	AUML
-----------	--	-----	------

**Auto Rerun SW**  ON  OFF

**Auto Rerun Range (result)**  ON  OFF

**Auto rerun condition (Absorbance)**

ABS. range	<input type="radio"/> ON <input type="radio"/> OFF
lower	<input type="radio"/> ON <input type="radio"/> OFF
high	<input type="radio"/> ON <input type="radio"/> OFF
Prozone range	<input type="radio"/> ON <input type="radio"/> OFF

	lower	higher
Serum	2	25
Urine		
Plasma		
CSF		
Dialysis		
other		

XXX: enter data by the user