


CHLORIDE

Ref. : CHLO-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	CHLO
DATA INFORMATION			
UNITS	mEq/L		
DECIMALS	1		
ANALYSIS			
TYPE	END		
Main W.Length 1	505		
Sub W.Length 1			
METHOD			
CORR		SLOPE	INTER
Y=	1	X+	0
CALIBRATION			
TYPE	Linear		
#1	xxx	#4	
#2		#5	
#3		#6	
NORMAL RANGE			
	MALE		FEMALE
	Low	High	Low High
Serum	XXX	XXX	XXX XXX
Urine			
Plasma			
CSF			
Dialysis			
other			

Item Name		XXX	CHLO
ASPIRATION			
KIND	<input type="radio"/> single <input type="radio"/> double		
VOLUME			
SAMPLE	3	μL	
REAGENT 1 vol	350		
DATA PROCESS		READ	ABSORBANCE LIMIT
Start	End	Low	High
Main 25	27	-3.000	
Sub		High	2.000
FACTOR		Endpoint limit	2
Blank correction	1	Linear check (%)	0
Third mix	<input type="radio"/> OFF <input type="radio"/> ON		
R1 Blank	<input type="radio"/> Water-blank <input type="radio"/> R1-blank-1		
Dilution			
Diluent	<input type="radio"/> 99: Dil 1 <input type="radio"/> 100: Dil 2		
PROZONE CHECK			
	start	End	Limit %
First			low <input type="radio"/> high <input type="radio"/>
Second			<input type="radio"/> <input type="radio"/>
third			<input type="radio"/> <input type="radio"/>

Item name		XXX	CHLO
Auto Rerun SW			
<input type="radio"/> ON <input type="radio"/> OFF			
Auto Rerun Range (result)			
<input type="radio"/> ON <input type="radio"/> OFF			
		Auto rerun condition (Absorbance)	
		ABS. range	
		lower	<input type="radio"/> ON <input type="radio"/> OFF
		high	<input type="radio"/> ON <input type="radio"/> OFF
		Prozone range	
		<input type="radio"/> ON	<input type="radio"/> OFF
	lower	higher	
Serum	12	130	
Urine	12	130	
Plasma			
CSF			
Dialysis			
other			

XXX: enter data by the user