


CREATININE JAFFE

“IMPROVED” TWO REGENT PROCEDURE

Ref. : CRCO-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	CRCO
DATA INFORMATION			
UNITS	mg/dL		
DECIMALS	2		
ANALYSIS			
TYPE	RATE		
Main W.Length 1	505		
Sub W.Length 2			
METHOD			
CORR	SLOPE	INTER	
Y=	1	X+	-0.20
CALIBRATION			
TYPE	Linear		
Standard			
#1	xxx	#4	
#2		#5	
#3		#6	
NORMAL RANGE			
	MALE		FEMALE
	Low	High	Low High
Serum	XXX	XXX	XXX XXX
Urine			
Plasma			
CSF			
Dialysis			
other			
Item Name		XXX	TIBC
ASPIRATION			
KIND	<input type="radio"/> single <input checked="" type="radio"/> double		
VOLUME			
Sample	30		
REAGENT 1 vol	150		
REAGENT 2 vol	150		
μL			
DATA PROCESS			
READ			
	Start	End	
Main	36	50	
Sub			
Absorbance limit			
Low	-3.000		
High	1.000		
FACTOR			
Blank correction			
Endpoint limit			
2			
Linear check (%)			
0			
Third mix	<input type="radio"/> OFF <input type="radio"/> ON		
R1 Blank	<input type="radio"/> Water-blank <input type="radio"/> R1-blank-1		
Dilution			
Diluent	<input type="radio"/> 99:Di 1 <input type="radio"/> 100:Di 2		
PROZONE CHECK			
	start	End	Limit %
First			
Second			
third			
			low high
			<input type="radio"/> <input type="radio"/>
			<input type="radio"/> <input type="radio"/>
Item name		XXX	CRCO
Auto Rerun SW			
<input type="radio"/> ON <input type="radio"/> OFF			
Auto Rerun Range (result)			
<input type="radio"/> ON <input type="radio"/> OFF			
Auto rerun condition (Absorbance)			
ABS. range			
lower <input type="radio"/> ON <input type="radio"/> OFF			
high <input type="radio"/> ON <input type="radio"/> OFF			
Prozone range <input type="radio"/> ON <input type="radio"/> OFF			
	lower	higher	
Serum	xxx	xxx	
Urine			
Plasma			
CSF			
Dialysis			
other			

XXX: enter data by the user