

# GLUCOSE HK SL

## TWO REAGENT PROCEDURE

Ref. : GHSL-XXXX

### APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

**Instruction for use** : 

**Working temperature** : 37°C

Item Name		XXX	GHSL				
<b>DATA INFORMATION</b>							
UNITS	mg/dL						
DECIMALS	0						
<b>ANALYSIS</b>							
TYPE	END						
Main W.Length 1	340						
Sub W.Length 2							
METHOD							
CORR	SLOPE	INTER					
Y=	1	X+	0				
<b>CALIBRATION</b>							
TYPE	Linear						
Standard							
#1	xxx	#4					
#2		#5					
#3		#6					
<b>NORMAL RANGE</b>							
	MALE		FEMALE				
	Low	High	Low	High			
Serum	XXX	XXX	XXX	XXX			
Urine							
Plasma							
CSF							
Dialysis							
other							
Item Name				XXX	GHSL		
<b>ASPIRATION</b>							
KIND	<input type="radio"/> single <input checked="" type="radio"/> double						
Sample	VOLUME						
REAGENT 1 vol	280	μL					
REAGENT 2 vol	70						
<b>DATA PROCESS</b>							
READ							
Main	Start	End	Absorbance limit				
Sub	52	54	Low	-3.000			
FACTOR			High	3.000			
Blank correction	1		Endpoint limit	2			
Dilution			Linear check (%)	0			
Diluent	<input checked="" type="radio"/> 99:Dil 1 <input type="radio"/> 100:Dil 2						
<b>PROZONE CHECK</b>							
Monitor		start	End	Limit %	low    high		
0 level point	1				<input type="radio"/> <input type="radio"/>		
SPAN	3				<input type="radio"/> <input type="radio"/>		
First					<input type="radio"/> <input type="radio"/>		
Second					<input type="radio"/> <input type="radio"/>		
third					<input type="radio"/> <input type="radio"/>		
Item name						XXX	GHSL
<b>Auto Rerun SW</b>							
<input checked="" type="radio"/> ON <input type="radio"/> OFF							
<b>Auto Rerun Range (result)</b>							
<input checked="" type="radio"/> ON <input type="radio"/> OFF <input type="radio"/> ON <input type="radio"/> OFF							
Auto rerun condition (Absorbance)		ABS. range					
		lower	<input checked="" type="radio"/> ON	<input type="radio"/> OFF			
		high	<input checked="" type="radio"/> ON	<input type="radio"/> OFF			
		Prozone range	<input type="radio"/> ON	<input type="radio"/> OFF			
Serum	xxx	xxx					
Urine							
Plasma							
CSF							
Dialysis							
other							

XXX: enter data by the user