


GLUCOSE PAP SL

Réf.: GPSL-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	GPSL
-----------	--	-----	------

DATA INFORMATION		CALIBRATION	
UNITS	mg/dL	TYPE	Linear
DECIMALS	0	Standard	
ANALYSIS		#1	xxx #4
TYPE	END	#2	#5
Main W.Length 1	505	#3	#6
Sub W.Length 2			
METHOD			
CORR		NORMAL RANGE	
SLOPE	1		
INTER	0		
Y=	X+		

Item Name		XXX	GPSL
-----------	--	-----	------

ASPIRATION		DATA PROCESS	
KIND	<input type="radio"/> single <input type="radio"/> double	READ	
VOLUME	3 μ L	Start	End
SAMPLE	3	Main	52 54
REAGENT 1 vol	300	Sub	
Third mix		Absorbance limit	
<input type="radio"/> OFF <input type="radio"/> ON		Low	-3.000
R1 Blank		High	2.000
<input type="radio"/> Water-blank <input type="radio"/> R1-blank-1		Endpoint limit 2	
Dilution		Linear check (%) 0	
Diluent			
<input type="radio"/> 99:Di1 <input type="radio"/> 100:Di2			
PROZONE CHECK			
Monitor		start	End
0 level point	1	Limit %	
SPAN	3	low	high
		<input type="radio"/>	<input type="radio"/>
		<input type="radio"/>	<input type="radio"/>

Item name		XXX	GPSL
-----------	--	-----	------

Auto Rerun SW		Auto rerun condition (Absorbance)	
<input type="radio"/> ON <input type="radio"/> OFF		ABS. range	
Auto Rerun Range (result)		lower	<input type="radio"/> ON <input type="radio"/> OFF
<input type="radio"/> ON <input type="radio"/> OFF		high	<input type="radio"/> ON <input type="radio"/> OFF
		Prozone range	<input type="radio"/> ON <input type="radio"/> OFF

XXX: enter data by the user