


IgA IP

Ref: IIGA-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

PROGRAM a Predilution 1/5
Calibrator/controls /samples

Item Name		XXX	IgA
DATA INFORMATION			
UNITS	mg/dL		
DECIMALS	2		
ANALYSIS			
TYPE	END		
Main W.Length 1	340		
Sub W.Length 2	800		
METHOD			
CORR		SLOPE	INTER
Y=	1	X+	0

CALIBRATION			
TYPE		Exponential	
#1	0.00	#4	xxx
#2	xxx	#5	xxx
#3	xxx	#6	xxx

NORMAL RANGE			
	MALE		FEMALE
	Low	High	High
Serum	xxx	xxx	xxx
Urine			
Plasma			
CSF			
Dialysis			
other			

Item Name		XXX	IgA
ASPIRATION			
KIND	<input type="radio"/> single <input checked="" type="radio"/> double		
SAMPLE	VOLUME		
REAGENT 1	5.6	μL	
vol	200		
REAGENT 2	33.6		
vol			
DATA PROCESS			
READ	Start	End	
Main	51	53	
Sub	29	31	
ABSORBANCE LIMIT			
	Low	-3.000	
	High	3.000	
FACTOR			
Blank correction	0.8595		
Endpoint limit		2	
Linear check (%)			
Third mix <input type="radio"/> OFF <input type="radio"/> ON			
R1 Blank <input checked="" type="radio"/> Water-blank <input type="radio"/> R1-blank-1			
Dilution			
Diluent	<input checked="" type="radio"/> 99:DI 1 <input type="radio"/> 100:DI 2		
PROZONE CHECK			
	start	End	Limit %
First			low <input type="radio"/> high <input type="radio"/>
Second			<input type="radio"/>
third			<input type="radio"/>

Item name	XXX	IgA							
Auto rerun condition (Absorbance)									
Auto Rerun SW	<input type="radio"/> ON <input type="radio"/> OFF		ABS. range		<input type="radio"/> ON	<input type="radio"/> OFF			
			lower						
			high						
Auto Rerun Range (result)									
	<input type="radio"/> ON	<input type="radio"/> OFF	<input type="radio"/> ON	<input type="radio"/> OFF	Prozone range		<input type="radio"/> ON	<input type="radio"/> OFF	
			lower	higher					
Serum	xxx	xxx							
Urine									
Plasma									
CSF									
Dialysis									
other									

XXX: enter data by the user

XXX: enter data by the user