


IgG IP

Ref: IIGG-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

PROGRAM a Predilution 1/20
Calibrator/controls /samples

Item Name XXX IgG

DATA INFORMATION

UNITS mg/dL
DECIMALS 2

ANALYSIS

TYPE END

Main W.Length 1 340
Sub W.Length 2 800

METHOD

CORR SLOPE INTER
Y= 1 X+ 0

CALIBRATION

TYPE Exponential

#1	0.00	#4	xxx
#2	xxx	#5	xxx
#3	xxx	#6	xxx

NORMAL RANGE

	MALE		FEMALE	
	Low	High	Low	High
Serum	xxx	xxx	xxx	xxx
Urine				
Plasma				
CSF				
Dialysis				
other				

Item Name XXX IgG

ASPIRATION

KIND single double

SAMPLE	VOLUME	
REAGENT 1	200	µL
REAGENT 2	20	

Third mix OFF ON

R1 Blank Water-blank R1-blank-1

Monitor

0 level point	1
SPAN	3

DATA PROCESS

READ	Start	End
Main	51	53
Sub	29	31

ABSORBANCE LIMIT

Low	-3.000
High	3.000

FACTOR

Blank correction 0.9104

Endpoint limit 2

Linear check (%)

Dilution

Diluent 99:DI 1 100:DI 2

PROZONE CHECK

	start	End	Limit %	low	high
First				<input type="radio"/>	<input type="radio"/>
Second				<input type="radio"/>	<input type="radio"/>
third				<input type="radio"/>	<input type="radio"/>

Item name	XXX	IgG							
Auto Rerun SW	<input type="radio"/> ON <input type="radio"/> OFF								
Auto Rerun Range (result)	<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF						
	lower	higher							
Serum	xxx	xxx							
Urine									
Plasma									
CSF									
Dialysis									
other									

XXX: enter data by the user

XXX: enter data by the user