


# IgM IP

Ref: IIGM-XXXX

**APPLICATION PRESTIGE 24i - PROPOSAL**

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

**PROGRAM a Predilution 1/5**  
**Calibrator/controls /samples**

---

Item Name: XXX IgM

**DATA INFORMATION**

UNITS: mg/dL  
DECIMALS: 2

**ANALYSIS**

TYPE: END  
Main W.Length 1: 340  
Sub W.Length 2: 800

METHOD:

CORR: SLOPE: INTER  
Y= 1 X+ 0

**CALIBRATION**

TYPE: Exponential

#1	0.00	#4	xxx
#2	xxx	#5	xxx
#3	xxx	#6	xxx

**NORMAL RANGE**

	MALE		FEMALE	
	Low	High	Low	High
Serum	xxx	xxx	xxx	xxx
Urine				
Plasma				
CSF				
Dialysis				
other				

---

Item Name: XXX IgM

**ASPIRATION**

KIND:  single  double

SAMPLE	VOLUME	μL
REAGENT 1 vol	200	
REAGENT 2 vol	24	

Third mix:  OFF  ON

R1 Blank:  Water-blank  R1-blank-1

Monitor: 0 level point: 1 SPAN: 3

**DATA PROCESS**

READ	Start	End
Main	51	53
Sub	29	31

ABSORBANCE LIMIT: Low: -3.000 High: 3.000

FACTOR: Blank correction: 0.8956 Endpoint limit: 2 Linear check (%):

Dilution: Diluent:  99:DI 1  100:DI 2

**PROZONE CHECK**

	start	End	Limit %	low	high
First				<input type="radio"/>	<input type="radio"/>
Second				<input type="radio"/>	<input type="radio"/>
third				<input type="radio"/>	<input type="radio"/>

---

Item name	XXX	IgM	Auto rerun condition (Absorbance)			
Auto Rerun SW			ABS. range			
<input checked="" type="radio"/> ON	<input type="radio"/> OFF		lower	<input checked="" type="radio"/> ON	<input type="radio"/> OFF	
			high	<input type="radio"/> ON	<input type="radio"/> OFF	
Auto Rerun Range (result)			Prozone range			
<input checked="" type="radio"/> ON	<input type="radio"/> OFF	<input type="radio"/> ON	<input type="radio"/> OFF	<input type="radio"/> ON	<input type="radio"/> OFF	
			lower	higher		
Serum	xxx	xxx				
Urine						
Plasma						
CSF						
Dialysis						
other						

XXX: enter data by the user

XXX: enter data by the user