


μALBUMIN IP

REF: IMAL-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

For more details, see the instruction for use : 

Working temperature : 37°C

PROCEDURE

Preparation of working reagent :

R1 : Reagent 1 Ref :IMAL-6125 Ready to use
R2 : Reagent 2 Ref :IMAL-5025 Ready to use

☞ Calibration :

μALBUMIN IP CALIBRATOR SET Ref : IMAL-0043
Ready to use.
Add a zero point (NaCl9g/L).

Item Name		XXX	IMAL	
DATA INFORMATION				
UNITS	mg/dL			
DECIMALS	2			
ANALYSIS				
TYPE	END			
Main W. Length 1	340			
Sub W. Length 2	700			
METHOD				
CORR	SLOPE	INTER		
Y=	1	X+	0	
CALIBRATION				
TYPE Exponential				
#1	0.00	#4	xxx	
#2	xxx	#5	xxx	
#3	xxx	#6	xxx	
NORMAL RANGE				
	MALE		FEMALE	
	Low	High	Low	High
Serum	xxx	xxx	xxx	xxx
Urine				
Plasma				
CSF				
Dialysis				
other				

Item Name		XXX	IMAL
ASPIRATION			
KIND	<input type="radio"/> single <input checked="" type="radio"/> double		
SAMPLE	VOLUME	μL	
REAGENT 1 vol	250		
REAGENT 2 vol	40		
DATA PROCESS			
READ	Start	End	
Main	52	54	
Sub	29	31	
ABSORBANCE LIMIT			
Low	-3.000		
High	3.000		
FACTOR			
Blank correction	0.8693		
Endpoint limit		2	
Linear check (%)			
Third mix	<input type="radio"/> OFF <input type="radio"/> ON		
R1 Blank	<input type="radio"/> Water-blank <input type="radio"/> R1-blank-1		
Dilution			
Diluent	<input checked="" type="radio"/> 99-Dil 1 <input type="radio"/> 100-Dil 2		
PROZONE CHECK			
	start	End	Limit %
First			
Second			
third			
			low <input type="radio"/> high <input type="radio"/>
			<input type="radio"/> <input type="radio"/>

Item name	XXX	IMAL							
Auto rerun condition (Absorbance)									
Auto Rerun SW	<input checked="" type="radio"/> ON	<input type="radio"/> OFF							
ABS. range									
								lower <input type="radio"/> higher <input type="radio"/>	
								<input type="radio"/> ON <input type="radio"/> OFF	
Auto Rerun Range (result)									
	<input checked="" type="radio"/> ON	<input type="radio"/> OFF							
								Prozone range <input type="radio"/> ON <input type="radio"/> OFF	
Serum	xxx	xxx							
Urine									
Plasma									
CSF									
Dialysis									
other									

XXX: enter data by the user

☞ Important modification from the previous version