

# RHEUMATOID FACTOR

REF: IRFA-XXXX

## APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use :  
 Consult Instruction for Use

<p>Item Name      RF</p> <p><u>DATA INFORMATION</u></p> <p>Units            IU/mL Decimals        0</p> <p><u>ANALYSIS</u></p> <p>Type            END</p> <p>W.Length 1      660 W.Length 2</p> <p>Method          Turbilatex</p> <p><u>CORR</u></p> <p>SLOPE           INTER 1.000 x +        0</p>	<p><u>CALIBRATION</u></p> <p>TYPE            Log2</p> <p>STANDARD</p> <table style="width: 100%; border: none;"> <tr> <td>#1</td><td>1/16 Cal. Val.</td><td>#4</td><td>1/2 Cal. Val.</td> </tr> <tr> <td>#2</td><td>1/8 Cal. Val.</td><td>#5</td><td>Cal. Val.</td> </tr> <tr> <td>#3</td><td>1/4 Cal. Val..</td><td>#6</td><td></td> </tr> </table> <p><u>NORMAL RANGE</u></p> <p>SERUM    MALE    LOW    HIGH              FEMALE</p> <p>URINE</p>	#1	1/16 Cal. Val.	#4	1/2 Cal. Val.	#2	1/8 Cal. Val.	#5	Cal. Val.	#3	1/4 Cal. Val..	#6																																	
#1	1/16 Cal. Val.	#4	1/2 Cal. Val.																																										
#2	1/8 Cal. Val.	#5	Cal. Val.																																										
#3	1/4 Cal. Val..	#6																																											
<p>Item Name      RF</p> <p><u>ASPIRATION</u></p> <p>KIND            Single <input type="checkbox"/> Double</p> <p>VOLUME**</p> <table style="width: 100%; border: none;"> <tr> <td>SAMPLE</td><td>3</td><td>µL</td> </tr> <tr> <td>REAGENT 1</td><td>240</td><td>µL</td> </tr> <tr> <td>REAGENT 2</td><td>60</td><td>µL</td> </tr> </table> <p>Third Mix.      <input checked="" type="checkbox"/> OFF   ON</p> <p>R1 Blank        Water   <input checked="" type="checkbox"/> R1-B</p> <p><u>MONITOR</u></p> <p>0 LEVEL POINT    1 SPAN                3.000</p>	SAMPLE	3	µL	REAGENT 1	240	µL	REAGENT 2	60	µL	<p><u>DATA PROCESS</u></p> <p><u>READ</u></p> <table style="width: 100%; border: none;"> <tr> <td></td><td>START</td><td>END</td><td>LOW</td><td>-3.000</td> </tr> <tr> <td>MAIN</td><td>41</td><td>42</td><td>HIGH</td><td>3.000</td> </tr> <tr> <td>SUB</td><td>33</td><td>34</td><td></td><td></td> </tr> </table> <p>ENDPOINT LIMIT    3</p> <p>LINEAR CHECK (%)</p> <p><u>FACTOR</u></p> <p>Blank Correction    1.000</p> <p><u>PROZONE CHECK</u></p> <table style="width: 100%; border: none;"> <tr> <td></td><td>START</td><td>END</td><td>LIMIT (%)</td><td></td> </tr> <tr> <td>FIRST</td><td></td><td></td><td></td><td></td> </tr> <tr> <td>SECOND</td><td></td><td></td><td></td><td><input checked="" type="checkbox"/> Low High</td> </tr> <tr> <td>THIRD</td><td></td><td></td><td></td><td><input checked="" type="checkbox"/> Low High</td> </tr> </table>		START	END	LOW	-3.000	MAIN	41	42	HIGH	3.000	SUB	33	34				START	END	LIMIT (%)		FIRST					SECOND				<input checked="" type="checkbox"/> Low High	THIRD				<input checked="" type="checkbox"/> Low High
SAMPLE	3	µL																																											
REAGENT 1	240	µL																																											
REAGENT 2	60	µL																																											
	START	END	LOW	-3.000																																									
MAIN	41	42	HIGH	3.000																																									
SUB	33	34																																											
	START	END	LIMIT (%)																																										
FIRST																																													
SECOND				<input checked="" type="checkbox"/> Low High																																									
THIRD				<input checked="" type="checkbox"/> Low High																																									

\*\* Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.