


# CHOLESTEROL LDL SL 2G

Ref. : LDL-XXXX

## APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	LDL
-----------	--	-----	-----

DATA INFORMATION		CALIBRATION	
UNITS	mg/dL	TYPE	Linear
DECIMALS	1	Standard	
ANALYSIS		#1	xxx #4
TYPE	END	#2	#5
Main W.Length 1	546	#3	#6
Sub W.Length 2		NORMAL RANGE	
METHOD			
CORR			
SLOPE			
Y-	1		
INTER			
X+	0		

Item Name		XXX	LDL
-----------	--	-----	-----

ASPIRATION		DATA PROCESS	
KIND	<input type="radio"/> single <input checked="" type="radio"/> double	READ	Absorbance limit
VOLUME	3	Start	End
Sample	295	Main	52 54
REAGENT 1 vol	100	Sub	29 31
REAGENT 2 vol		FACTOR	Blank correction
Third mix			0.74874
<input type="radio"/> OFF <input type="radio"/> ON		Endpoint limit	2
R1 Blank	<input type="radio"/> Water-blank <input type="radio"/> R1-blank-1	Linear check (%)	0
Dilution		PROZONE CHECK	
Diluent	<input type="radio"/> 99: Dil 1 <input type="radio"/> 100: Dil 2	start	End
Monitor		Limit %	
0 level point	1	First	
SPAN	3	Second	
		third	
		low	high
		<input type="radio"/>	<input type="radio"/>
		<input type="radio"/>	<input type="radio"/>

Item name		XXX	LDL
-----------	--	-----	-----

Auto Rerun SW		Auto rerun condition (Absorbance)	
<input type="radio"/> ON <input type="radio"/> OFF		ABS. range	
Auto Rerun Range (result)		lower	<input type="radio"/> ON <input type="radio"/> OFF
<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF	high	<input type="radio"/> ON <input type="radio"/> OFF
		Prozone range	<input type="radio"/> ON <input type="radio"/> OFF
lower	higher		
Serum	15 700		
Urine			
Plasma			
CSF			
Dialysis			
other			

XXX: enter data by the user