


ALP (DEA) SL

TWO REAGENT PROCEDURE

Ref.: PASL-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	PASL
-----------	--	-----	------

DATA INFORMATION		CALIBRATION	
UNITS	U/L	TYPE	Factor
DECIMALS	0	Factor	3424

ANALYSIS		NORMAL RANGE	
TYPE	RATE		
Main W.Length 1	405		
Sub W.Length 2			
METHOD	DEA		

CORR	SLOPE	INTER
Y=	1	X+
		0

Item Name		XXX	PASL
-----------	--	-----	------

ASPIRATION		DATA PROCESS		ABSORBANCE LIMIT	
KIND	<input type="radio"/> single <input checked="" type="radio"/> double	READ	Start	End	
SAMPLE	VOLUME	Main	39	44	Low
REAGENT 1 vol	4	Sub			High
REAGENT 2 vol	200				
	50				
	µL				
FACTOR		Endpoint limit		0	
Blank correction	1	Linear check		40	
		Linear check (%)			

Third mix	<input type="radio"/> OFF <input type="radio"/> ON				
R1 Blank	<input type="radio"/> Water-blank <input type="radio"/> R1-blank-1				
Dilution					
Diluent	<input type="radio"/> 99-Dil 1 <input type="radio"/> 100-Dil 2				
PROZONE CHECK					
	start	End	Limit %	low	high
First				<input type="radio"/>	<input type="radio"/>
Second				<input type="radio"/>	<input type="radio"/>
third				<input type="radio"/>	<input type="radio"/>

Item name		XXX	PASL
-----------	--	-----	------

Auto Rerun SW		Auto rerun condition (Absorbance)	
<input checked="" type="radio"/> ON <input type="radio"/> OFF		ABS. range	<input type="radio"/> lower <input type="radio"/> ON <input type="radio"/> OFF
		high	<input type="radio"/> ON <input type="radio"/> OFF
Auto Rerun Range (result)		Prozone range	
<input type="radio"/> ON <input type="radio"/> OFF	<input type="radio"/> ON <input type="radio"/> OFF		<input type="radio"/> ON <input type="radio"/> OFF

	lower	higher
Serum	30	900
Urine		
Plasma		
CSF		
Dialysis		
other		

XXX: enter data by the user