


PHOSPHORUS

Réf. : PHOS-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

Item Name		XXX	PHOS
DATA INFORMATION			
UNITS	mg/dL		
DECIMALS	2		
ANALYSIS			
TYPE	END		
Main W.Length 1	340		
Sub W.Length 1			
METHOD			
CORR	SLOPE	INTER	
Y=	1	X+	0
CALIBRATION			
TYPE	Linear		
#1	xxx	#4	
#2		#5	
#3		#6	
NORMAL RANGE			
	MALE		FEMALE
	Low	High	Low High
Serum	XXX	XXX	XXX XXX
Urine			
Plasma			
CSF			
Dialysis			
other			

Item Name		XXX	PHOS
ASPIRATION			
KIND	<input type="radio"/> single <input type="radio"/> double		
VOLUME			
SAMPLE	3	µL	
REAGENT 1 vol	300		
DATA PROCESS			
READ	Start	End	
Main	35	37	
Sub			
ABSORBANCE LIMIT			
Low	-3.000		
High	2.000		
FACTOR			
Blank correction	1		
Endpoint limit			
Endpoint limit	2		
Linear check (%)	0		
Third mix <input type="radio"/> OFF <input type="radio"/> ON			
R1 Blank <input type="radio"/> Water-blank <input type="radio"/> R1-blank-1			
Dilution			
Diluent	<input type="radio"/> 99:Di1 <input type="radio"/> 100:Di2		
PROZONE CHECK			
	start	End	Limit %
First			low high
Second			<input type="radio"/> <input type="radio"/>
third			<input type="radio"/> <input type="radio"/>

Item name		XXX	PHOS
Auto Rerun SW			
<input type="radio"/> ON <input type="radio"/> OFF			
Auto Rerun Range (result)			
<input type="radio"/> ON <input type="radio"/> OFF			
Auto rerun condition (Absorbance)			
ABS range		<input type="radio"/> ON <input type="radio"/> OFF	
lower		<input type="radio"/> ON	<input type="radio"/> OFF
high		<input type="radio"/> ON	<input type="radio"/> OFF
Prozone range		<input type="radio"/> ON <input type="radio"/> OFF	
Monitor			
0 level point	1		
SPAN	3		
lower higher			
Serum	2	20	
Urine			
Plasma			
CSF			
Dialysis			
other			

XXX: enter data by the user