


TRIGLYCERIDES MONO SL NEW

Réf.: TGML-XXXX

APPLICATION PRESTIGE 24i - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use

Instruction for use : 

Working temperature : 37°C

| | | | |
|-----------|--|-----|------|
| Item Name | | XXX | TGML |
|-----------|--|-----|------|

DATA INFORMATION

| | |
|----------|-------|
| UNITS | mg/dL |
| DECIMALS | 0 |

ANALYSIS

| | |
|-----------------|-----|
| TYPE | END |
| Main W.Length 1 | 505 |
| Sub W.Length 1 | |

METHOD

| | | |
|------|-------|-------|
| CORR | SLOPE | INTER |
| Y= | 1 | X= 0 |

CALIBRATION

| | |
|----------|--------|
| TYPE | Linear |
| Standard | |
| #1 | xxx #4 |
| #2 | #5 |
| #3 | #6 |

NORMAL RANGE

| | MALE | | FEMALE | |
|----------|------|------|--------|------|
| | Low | High | Low | High |
| Serum | XXX | XXX | XXX | XXX |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| other | | | | |

| | | | |
|-----------|--|-----|------|
| Item Name | | XXX | TGML |
|-----------|--|-----|------|

Aspiration

KIND single double

| | |
|----------------------|--------|
| VOLUME | |
| SAMPLE REAGENT 1 vol | 300 µL |

Data Process

Read

| | | |
|------|-------|-----|
| Main | Start | End |
| Sub | 42 | 44 |

Absorbance limit

| | |
|------|--------|
| Low | -3.000 |
| High | 1.800 |

Blank correction 1

Endpoint limit 2

Linear check (%) 0

Third mix OFF ON

R1 Blank Water-blank R1-blank-1

Dilution

Diluent 99-Dil 1 100-Dil 2

PROZONE CHECK

| | start | End | Limit % | low | high |
|--------|-------|-----|---------|-----------------------|-----------------------|
| First | | | | <input type="radio"/> | <input type="radio"/> |
| Second | | | | <input type="radio"/> | <input type="radio"/> |
| third | | | | <input type="radio"/> | <input type="radio"/> |

| | | | |
|-----------|--|-----|------|
| Item name | | XXX | TGML |
|-----------|--|-----|------|

Auto Rerun SW

ON OFF

Auto Rerun Range (result)

ON OFF

Auto rerun condition (Absorbance)

ABS. range lower ON OFF

high ON OFF

Prozone range ON OFF

| | lower | higher |
|----------|-------|--------|
| Serum | 20 | 1000 |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| other | | |

XXX: enter data by the user