

## APPLICATION - RESPONS 920 - ERBA XL200 / XL-180 - EM200 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Consult Instruction for Use

### TEST DETAILS

Test	:	ALBU	Auto Rerun	<input type="checkbox"/>
Report Name	:	ALBUMIN	Online calibration	<input type="checkbox"/>
Unit	:	g/dL	Cuvette Wash	<input type="checkbox"/>
Wavelength-Primary	:	600	Total reagents	: 1
Assay Type	:	1-Point	Curve Type	: Linear
M1 Start	:	0	M1 End	: 0
M2 Start	:	3	M2 End	: 3
Sample Replicates	:	1	Standard replicates	: 3
Control Replicates	:	1	Control Interval	: 0
Reaction Direction	:	Increasing	React. Abs. Limit	: *
Prozone Limit %	:	0	Prozone Check	: lower
Linearity Limit %	:	0	Delta Abs./Min.	: 0.0000
Technical Minimum	:	*	Technical Maximum	: *
Y= aX + b	a=	: 1.0000	B=	: 0.0000

### TEST VOLUMES

Test : ALBU  
Sample Type : Serum

#### Sample Volumes

Normal	:	2.00 µL	Dilution Ratio	:	1 X
Increase	:	* µL	Dilution Ratio	:	* X
Decrease	:	* µL	Dilution Ratio	:	* X
Standard Volume	:	2.00 µL			

#### reagent Volumes and Stirrer Speed

RGT-1 Volume	:	300 µL	R1 Stirrer Speed	:	high
RGT-2 Volume	:	0 µL	R2 Stirrer Speed	:	0

#### Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

### REFERENCE RANGES

Test : ALBU  
Sample Type : Serum  
Reference Range : DEFAULT  
Category : \*

#### Reference Range

		Lower Limit	Lower Limit
		g/dL	g/dL
Normal	:	*	*
Panic	:	*	*

#### Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

\*: Values entered by the operator