

APPLICATION - RESPONS 920 - ERBA XL200 / XL-180 - EM200 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Consult Instruction for Use

TEST DETAILS

Test	:	ASO			Auto Rerun	<input type="checkbox"/>		
Report Name	:	ANTI-STREPTOLYSIN O			Online calibration	<input type="checkbox"/>		
Unit	:	IU/mL	Decimal Places	:	1	Cuvette Wash	<input type="checkbox"/>	
Wavelength-Primary	:	546	Secondary	:	-	Total reagents	:	2
Assay Type	:	2-Point	Curve Type	:	Linear	Reagent 1	:	ASLO R1
M1 Start	:	14	M1 End	:	15	Reagent 2	:	ASLO R2
M2 Start	:	24	M2 End	:	25			
Sample Replicates	:	1	Standard replicates	:	3	Consumables/Calibrators:		
Control Replicates	:	1	Control Interval	:	0	Blank/Level 0	:	0
Reaction Direction	:	Increasing	React. Abs. Limit	:	*	Calibrator 1	:	*
Prozone Limit %	:	0	Prozone Check	:	lower			
Linearity Limit %	:	0	Delta Abs./Min.	:	0.0000			
Technical Minimum	:	*	Technical Maximum	:	*			
Y= aX + b	a=	1.0000	B=	:	0.0000			

TEST VOLUMES

Test : ASO
Sample Type : Serum

Sample Volumes

Normal	:	3 µL	Dilution Ratio	:	1 X
Increase	:	* µL	Dilution Ratio	:	* X
Decrease	:	* µL	Dilution Ratio	:	* X
Standard Volume	:	3 µL			

reagent Volumes and Stirrer Speed

RGT-1 Volume	:	240 µL	R1 Stirrer Speed	:	high
RGT-2 Volume	:	60 µL	R2 Stirrer Speed	:	high

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

REFERENCE RANGES

Test : ASO
Sample Type : Serum
Reference Range : DEFAULT
Category : *

Reference Range

		Lower Limit		Lower Limit
		IU/mL		IU/mL
Normal	:	*	:	*
Panic	:	*	:	*

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

*: Values entered by the operator