

APPLICATION - RESPONS 920 - ERBA XL200 / XL-180 - EM200 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Consult Instruction for Use

TEST DETAILS

Test	: CKMB	Auto Rerun	<input type="checkbox"/>
Report Name	: CK-MB SL	Online calibration	<input type="checkbox"/>
Unit	: U/L	Decimal Places	: 0
Wavelength-Primary	: 340	Secondary	: -
Assay Type	: RATE-A	Curve Type	: K-FACTOR
M1 Start	: 0	M1 End	: 0
M2 Start	: 12	M2 End	: 28
Sample Replicates	: 1	Factor	: 4133
Control Replicates	: 1	Control Interval	: 0
Reaction Direction	: Increasing	React. Abs. Limit	: *
Prozone Limit %	: 0	Prozone Check	: lower
Linearity Limit %	: 0	Delta Abs./Min.	: 0.0000
Technical Minimum	: *	Technical Maximum	: *
Y= aX + b	a= : 1.0000	B= :	0.0000
		Reagent 1	: CMSL WR
		Reagent 2	:

TEST VOLUMES

Test : CKMB
Sample Type : Serum

Sample Volumes

Normal	: 16.00 µL	Dilution Ratio	: 1 X
Increase	: * µL	Dilution Ratio	: * X
Decrease	: * µL	Dilution Ratio	: * X
Standard Volume	: 16.00 µL		

reagent Volumes and Stirrer Speed

RGT-1 Volume	: 200 µL	R1 Stirrer Speed	: medium
RGT-2 Volume	: 0 µL	R2 Stirrer Speed	: 0

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

REFERENCE RANGES

Test : CKMB
Sample Type : Serum
Reference Range : DEFAULT
Category : *

Reference Range

	Lower Limit	Lower Limit
	U/L	U/L
Normal	: *	: *
Panic	: *	: *

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

*: Values entered by the operator