

APPLICATION - RESPONS 920 - ERBA XL200 / XL-180 - EM200 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Consult Instruction for Use

TEST DETAILS

Test	: CRP	Auto Rerun	<input type="checkbox"/>
Report Name	: CRP IP	Online calibration	<input type="checkbox"/>
Unit	: g/dL	Decimal Places	: 1
Wavelength-Primary	: 340	Secondary	: -
Assay Type	: 2-Point	Curve Type	: Cubic spline
M1 Start	: 15	M1 End	: 15
M2 Start	: 31	M2 End	: 31
Sample Replicates	: 1	Standard replicates	: 3
Control Replicates	: 1	Control Interval	: 0
Reaction Direction	: Increasing	React. Abs. Limit	: *
Prozone Limit %	: 0	Prozone Check	: lower
Linearity Limit %	: 0	Delta Abs./Min.	: 0.0000
Technical Minimum	: *	Technical Maximum	: *
Y= aX + b	a= : 1.0000	B= :	0.0000
		Cuvette Wash	<input type="checkbox"/>
		Total reagents	: 2
		Reagent 1	: ICRP R1
		Reagent 2	: ICRP R2

TEST VOLUMES

Test : CRP
Sample Type : Serum

Sample Volumes

Normal	: 11.50 µL	Dilution Ratio	: 1 X
Increase	: * µL	Dilution Ratio	: * X
Decrease	: * µL	Dilution Ratio	: * X
Standard Volume	: 11.50 µL		

reagent Volumes and Stirrer Speed

RGT-1 Volume	: 180 µL	R1 Stirrer Speed	: high
RGT-2 Volume	: 18 µL	R2 Stirrer Speed	: high

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

REFERENCE RANGES

Test : CRP
Sample Type : Serum
Reference Range : DEFAULT
Category : *

Reference Range

	Lower Limit	Lower Limit
	g/dL	g/dL
Normal	: *	: *
Panic	: *	: *

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

*: Values entered by the operator