

## APPLICATION - RESPONS 920 - ERBA XL200 / XL-180 - EM200 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Consult Instruction for Use

### TEST DETAILS

Test	:	RF	Auto Rerun	<input type="checkbox"/>	
Report Name	:	RHEUMATOID FACTOR	Online calibration	<input type="checkbox"/>	
Unit	:	IU/mL	Cuvette Wash	<input type="checkbox"/>	
Wavelength-Primary	:	600	Secondary	:	-
Assay Type	:	2-Point	Curve Type	:	5P Logit-Log
M1 Start	:	15	M1 End	:	15
M2 Start	:	24	M2 End	:	24
Sample Replicates	:	1	Standard replicates	:	1
Control Replicates	:	1	Control Interval	:	0
Reaction Direction	:	Increasing	React. Abs. Limit	:	*
Prozone Limit %	:	0	Prozone Check	:	lower
Linearity Limit %	:	0	Delta Abs./Min.	:	0.0000
Technical Minimum	:	*	Technical Maximum	:	*
Y= aX + b	a=	1.0000	B=	:	0.0000
			Reagent 1	:	IRFA R1
			Reagent 2	:	IRFA R2
			Consumables/Calibrators:		
			Blank/Level 0	:	0
			Calibrator 1	:	*
			Calibrator 2	:	*
			Calibrator 3	:	*
			Calibrator 4	:	*
			Calibrator 5	:	*
			Calibrator 6	:	*

### TEST VOLUMES

Test : RF  
Sample Type : Serum

#### Sample Volumes

Normal	:	2.8 µL	Dilution Ratio	:	1 X
Increase	:	* µL	Dilution Ratio	:	* X
Decrease	:	* µL	Dilution Ratio	:	* X
Standard Volume	:	2.8 µL			

#### reagent Volumes and Stirrer Speed

RGT-1 Volume	:	220 µL	R1 Stirrer Speed	:	high
RGT-2 Volume	:	55 µL	R2 Stirrer Speed	:	high

#### Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

### REFERENCE RANGES

Test : RF  
Sample Type : Serum  
Reference Range : DEFAULT  
Category : \*

#### Reference Range

		Lower Limit	Lower Limit
		IU/mL	IU/mL
Normal	:	*	*
Panic	:	*	*

#### Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

\*: Values entered by the operator