

APPLICATION - RESPONS 920 - ERBA XL200 / XL-180 - EM200 - PROPOSAL

This application is intended to serve as a guide for using the referenced ELITechGroup Reagent on this instrument system only. It is recommended that the user validate this application prior to routine use.

Consult Instruction for Use

TEST DETAILS

Test	: PROT	Auto Rerun	<input type="checkbox"/>
Report Name	: TOTAL PROTEIN PLUS	Online calibration	<input type="checkbox"/>
Unit	: g/dL	Cuvette Wash	<input type="checkbox"/>
Wavelength-Primary	: 546	Decimal Places	: 1
Assay Type	: 1-Point	Secondary	: -
M1 Start	: 0	Curve Type	: Linear
M2 Start	: 35*	M1 End	: 0
Sample Replicates	: 1	M2 End	: 35*
Control Replicates	: 1	Standard replicates	: 3
Reaction Direction	: Increasing	Control Interval	: 0
Prozone Limit %	: 0	React. Abs. Limit	: *
Linearity Limit %	: 0	Prozone Check	: lower
Technical Minimum	: *	Delta Abs./Min.	: 0.0000
Y= aX + b	a= : 1.0000	Technical Maximum	: *
		B=	: 0.0000
		Reagent 1	: PROB R1
		Reagent 2	:

TEST VOLUMES

Test : PROT
Sample Type : Serum

Sample Volumes

Normal	: 2.00 µL	Dilution Ratio	: 1 X
Increase	: * µL	Dilution Ratio	: * X
Decrease	: * µL	Dilution Ratio	: * X
Standard Volume	: 2.00 µL		

reagent Volumes and Stirrer Speed

RGT-1 Volume	: 200 µL	R1 Stirrer Speed	: high
RGT-2 Volume	: 0 µL	R2 Stirrer Speed	: 0

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

REFERENCE RANGES

Test : PROT
Sample Type : Serum
Reference Range : DEFAULT
Category : *

Reference Range

	Lower Limit	Lower Limit
	g/dL	g/dL
Normal	: *	: *
Panic	: *	: *

Sample Types

- Serum
- Urine
- CSF
- Plasma
- Whole blood
- Other

* XL-180 : 34

*: Values entered by the operator