

A15 / A25

BioSYSTEMS

AA-A15/A25-CAT-1

A	Versions	
ACID PHOSPHATASE	PACI-1	PRO
ALBUMIN	ALBU-1	PRO
ALP (DEA) SL monoprocédure	PASLmono-1	PRO
ALP (DEA) SL biprocédure.....	PASLbi-1	PRO
ALT/GPT 4+1 SL monoprocédure.....	ALSL4+1mono-1..	PRO
ALT/GPT 4+1 SL biprocédure	ALSL4+1bi-1	PRO
AMYLASE SL.....	AMSL-1	PRO
AST/GOT 4+1 SL monoprocédure.....	ASSL4+1mono-1 .	PRO
AST/GOT 4+1 SL biprocédure.....	ASSL4+1bi-1	PRO

B		
BILIRUBIN TOTAL & DIRECT 4+1 ...	BIDT	PRO

C		
CALCIUM ARSENAZO.....	CALA	PRO
CALCIUM OCPC monoprocédure	CALOmono-1	PRO
CHLORIDE	CHLO-1	PRO
CHOLESTEROL SL	CHSL-1	PRO
CHOLESTEROL HDL	HDLC-1	PRO
CHOLESTEROL HDL SL 2G	HDLL-1	PRO
CHOLESTEROL LDL SL 2G	LDLL-1	PRO
CKMB SL monoprocédure	CMSLmono-1.....	PRO
CK NAC SL monoprocédure	CKSLmono-1	PRO
CREATININE JAFFE monoprocédure ..	CRCOmono-1	PRO
CREATININE JAFFE biprocédure	CRCObi-1	PRO

G		
GAMMA GT SL monoprocédure	GASLmono-1	PRO
GAMMA GT SL biprocédure.....	GASLbi-1	PRO
GLUCOSE HK SL monoprocédure ..	GHSLmono-1.....	PRO
GLUCOSE HK SL biprocédure.....	GHSLbi-1	PRO
GLUCOSE PAP SL	GPSL-1	PRO

I		
IRON FERROZINE	FEFR-1	PRO

L	Versions	
LACTATE.....	LACT-1	PRO
LDH-P 4+1 SL monoprocédure	LDSL4+1mono-1 .	PRO
LDH-P 4+1 SL biprocédure.....	LDSL4+1bi-1	PRO

M		
MAGNESIUM CALMAGITE	MAGNmono-1	PRO
MICROPROTEIN	PRTP-1	PRO


P		
PHOSPHORUS.....	PHOS-1	PRO

T		
TOTAL PROTEIN.....	PRTB-1	PRO
TRIGLYCERIDES MONO SL NEW....	TGMLN-1	PRO

U		
UREA UV SL monoprocédure	URSLmono-1	PRO
UREA UV SL biprocédure.....	URSLbi-1.....	PRO
URIC ACID MONO SL	AUML-1	PRO
URIC ACID SL monoprocédure.....	AUSLmono-1	PRO
URIC ACID SL biprocédure	AUSLbi-1	PRO

PROPOSAL (PRO): These applications have been established theoretically. It is advised to check them.


APPLICATION A15/A25**PROPOSAL**

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	PACI	PACI
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	25	25
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	405	405
Reference	-	-	
Times			
Reading 1 (sec/cycle)	288 sec / cycle n° 13	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	408 sec / cycle n° 18	405 sec / cycle n° 28	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	40	40	

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
Working temperature : 37°C

		A15	A25
GENERAL	Test name	ALBU	ALBU
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	g/dL	g/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	600	600
Reference	-	-	
Times			
Reading 1 (sec/cycle)	72 sec / cycle n° 4	75 sec / cycle n° 6	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	6	6	

xxx : User defined

REF: PASL-0400 2 x 62.5 mL
 PASL-0420 4 x 62.5 mL
 PASL-0500 5 x 125 mL

APPLICATION A15/A25

PROPOSAL

For more details, see the instructions for use: 


Working temperature : 37°C

		A15	A25
GENERAL	Test name	PASL	PASL
	Analysis mode	Monoreagent Kinetics	Monoreagent Kinetics
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	5	5
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	405	405
Reference	-	-	
Times	Reading 1 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5
	Reading 2 (sec/cycle)	192 sec / cycle n° 9	195 sec / cycle n° 14
	Reagent 2 (sec/cycle)	-	-
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	900	900	

xxx : User defined

REF: PASL-0400 2 x 62.5 mL
 PASL-0420 4 x 62.5 mL
 PASL-0500 5 x 125 mL

APPLICATION A15/A25
PROPOSAL


For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	PASL	PASL
	Analysis mode	Bireagent Kinetics	Bireagent Kinetics
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	5	5
	Reagent 1 (µL)	200	200
	Reagent 2 (µL)	50	50
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	405	405
Reference	-	-	
Times			
Reading 1 (sec/cycle)	96 sec / cycle n° 5	120 sec / cycle n° 9	
Reading 2 (sec/cycle)	240 sec / cycle n° 11	255sec / cycle n° 18	
Reagent 2 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	900	900	

xxx : User defined

REF: ALSL-0410 2 x 62.5 mL
 ALSL-0430 4 x 62.5 mL
 ALSL-0510 5 x 125 mL

APPLICATION A15/A25
PROPOSAL


For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	ALSL	ALSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	30	30
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
Reading 2 (sec/cycle)	216 sec / cycle n° 10	225 sec / cycle n° 16	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	250	250	

xxx : User defined

REF: ALSL-0410 2 x 62.5 mL
 ALSL-0430 4 x 62.5 mL
 ALSL-0510 5 x 125 mL


APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	ALSL	ALSL
	Analysis mode	Bireagent Kinetic	Bireagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report		xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	30	30
	Reagent 1 (µL)	240	240
	Reagent 2 (µL)	60	60
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	96 sec / cycle n° 5	120 sec / cycle n° 9	
Reading 2 (sec/cycle)	264 sec / cycle n° 12	285sec / cycle n° 20	
Reagent 2 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit		xxx	xxx
Kinetik Blank limit		-	-
Linearity Limit		250	250

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	AMSL	AMSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	405	405
Reference	-	-	
Times			
Reading 1 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
Reading 2 (sec/cycle)	216 sec / cycle n° 10	225 sec / cycle n° 16	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	2000	2000	

xxx : User defined



REF: ASSL-0410 2 x 62.5 mL
 ASSL-0430 4 x 62.5 mL
 ASSL-0510 5 x 125 mL

APPLICATION A15/A25	PROPOSAL
----------------------------	-----------------

For more details, see the instructions for use:
 Working temperature : 37°C


		A15	A25
GENERAL	Test name	ASSL	ASSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	30	30
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
Reading 2 (sec/cycle)	216 sec / cycle n° 10	225 sec / cycle n° 16	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	250	250

xxx : User defined

REF: ASSL-0410 2 x 62.5 mL
 ASSL-0430 4 x 62.5 mL
 ASSL-0510 5 x 125 mL

APPLICATION A15/A25

PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	ASSL	ASSL
	Analysis mode	Bireagent Kinetic	Bireagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	30	30
	Reagent 1 (µL)	240	240
	Reagent 2 (µL)	60	60
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	96 sec / cycle n° 5	120 sec / cycle n° 9	
Reading 2 (sec/cycle)	264 sec / cycle n° 12	285sec / cycle n° 20	
Reagent 2 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	250	250

xxx : User defined



BILIRUBIN TOTAL & DIRECT 4+1

REF: BITD-0600	Total & Direct 4+1	2 x 125 mL
BIDI- 0600	Direct 4+1	2 x 125 mL
BITO-0600	Total 4+1	2 x 125 mL

APPLICATION A15/A25

PROPOSAL


For more details, see the instructions for use:
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	BIDT	BIDT
	Analysis mode	Bireagent differential	Bireagent differential
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	30	30
	Reagent 1 (µL)	240	240
	Reagent 2 (µL)	60	60
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	560	560
Reference	-	-	
Times			
Reading 1 (sec/cycle)	264 sec / cycle n° 12	270 sec / cycle n° 19	
Reading 2 (sec/cycle)	600 sec / cycle n° 26	600 sec / cycle n° 41	
Reagent 2 (sec/cycle)	288 sec / cycle n° 13	285 sec / cycle n° 20	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	20 bili Total / 18 bili Direct	20 bili Total / 18 bili Direct	

xxx : User defined

(09/2009)
AA-A15/A25-BIDT-1


APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	CALA	CALA
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
Test name in patient report		xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	10	10
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	600	600
Reference	-	-	
Times			
Reading 1 (sec/cycle)	72 sec / cycle n° 4	75 sec / cycle n° 6	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit		xxx	xxx
Kinetik Blank limit		-	-
Linearity Limit		15	15

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	CALO	CALO
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	6	6
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	560	560
Reference	-	-	
Times			
Reading 1 (sec/cycle)	72sec / cycle n° 4	75 sec / cycle n° 6	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	13.5	13.5	

xxx : User defined

APPLICATION A15/A25
PROPOSAL

 For more details, see the instructions for use: 

Working temperature : 37°C


		A15	A25
GENERAL	Test name	CHLO	CHLO
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mEq/L	mEq/L
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
	Reference	-	-
Times			
Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	130	130	

xxx : User defined

REF:	
CHSL-0490	1 x 100 mL
CHSL-0500	6 X 100 mL
CHSL-0700	4 X 250 mL
CHSL-0507	6 X 100 mL +STD
CHSL-0707	4 X 250 mL +STD

APPLICATION A15/A25


PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	CHSL	CHSL
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
Reference	-	-	
Times			
Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	600	600	

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C


		A15	A25
GENERAL	Test name	HDLC	HDLC
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	15	15
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
Reference	-	-	
Times			
Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	200	200	

xxx : User defined

REF: HDLL-0380 3 x 80 mL
 HDLL-0390 3 x 80 mL

APPLICATION A15/A25


PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	HDLL	HDLL
	Analysis mode	Bireagent differential	Bireagent differential
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	100	100
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	600	600
Times	Reference	-	-
	Reading 1 (sec/cycle)	288 sec / cycle n° 13	285 sec / cycle n° 20
	Reading 2 (sec/cycle)	600 sec / cycle n° 26	600 sec / cycle n° 41
	Reagent 2 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	200	200	

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	LDLL	LDLL
	Analysis mode	Bireagent differential	Bireagent differential
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	100	100
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	560	560
Times	Reference	-	-
	Reading 1 (sec/cycle)	288 sec / cycle n° 13	285 sec / cycle n° 20
	Reading 2 (sec/cycle)	600 sec / cycle n° 26	600 sec / cycle n° 41
	Reagent 2 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	700	700	

xxx : User defined

APPLICATION A15/A25

PROPOSAL


For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	CMSL	CMSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	10	10
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	144 sec / cycle n° 7	135 sec / cycle n° 10	
Reading 2 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	600	600	

xxx : User defined

REF: CKSL-0410 2 x 62.5 mL
 CKSL-0430 4 x 62.5 mL


APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	CKSL	CKSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	10	10
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	144 sec / cycle n° 7	135 sec / cycle n° 10	
Reading 2 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	1700	1700	

xxx : User defined


APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
Working temperature : 37°C

		A15	A25
GENERAL	Test name	CRCO	CRCO
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	25	25
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
Times	Reference	-	-
	Reading 1 (sec/cycle)	24 sec / cycle n° 2	15 sec / cycle n° 2
	Reading 2 (sec/cycle)	144 sec / cycle n° 7	135 sec / cycle n° 10
	Reagent 2 (sec/cycle)	-	-
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	15	15	

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
Working temperature : 37°C


		A15	A25
GENERAL	Test name	CRCO	CRCO
	Analysis mode	Bireagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
Test name in patient report		xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	25	25
	Reagent 1 (µL)	125	125
	Reagent 2 (µL)	125	125
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
Reference	-	-	
Times	Reading 1 (sec/cycle)	48 sec / cycle n° 3	45 sec / cycle n° 4
	Reading 2 (sec/cycle)	168 sec / cycle n° 8	165 sec / cycle n° 12
	Reagent 2 (sec/cycle)	24 sec / cycle n° 2	15 sec / cycle n° 2
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit		xxx	xxx
Kinetik Blank limit		-	-
Linearity Limit		15	15

xxx : User defined

REF: GASL-0400 2 x 62.5 mL
 GASL-0420 4 x 62.5 mL
 GASL-0500 5 x 125 mL

APPLICATION A15/A25

PROPOSAL


For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	GASL	GASL
	Analysis mode	Monoreagent Kinetics	Monoreagent Kinetics
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	30	30
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	405	405
Reference	-	-	
Times			
Reading 1 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
Reading 2 (sec/cycle)	216 sec / cycle n° 10	225 sec / cycle n° 16	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	1000	1000	

xxx : User defined

REF: GASL-0400 2 x 62.5 mL
 GASL-0420 4 x 62.5 mL
 GASL-0500 5 x 125 mL


APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	GASL	GASL
	Analysis mode	Bireagent Kinetics	Bireagent Kinetics
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report		xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	30	30
	Reagent 1 (µL)	240	240
	Reagent 2 (µL)	60	60
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	405	405
Reference	-	-	
Times	Reading 1 (sec/cycle)	96 sec / cycle n° 5	120 sec / cycle n° 9
	Reading 2 (sec/cycle)	264 sec / cycle n° 12	285sec / cycle n° 20
	Reagent 2 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit		xxx	xxx
Kinetik Blank limit		-	-
Linearity Limit		1000	1000

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	GHSL	GHSL
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Times	Reference	-	-
	Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21
	Reading 2 (sec/cycle)	-	-
	Reagent 2 (sec/cycle)	-	-
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
	Kinetik Blank limit	-	
	Linearity Limit	600	

xxx : User defined

APPLICATION A15/A25
PROPOSAL

 For more details, see the instructions for use: 

Working temperature : 37°C

		A15	A25
GENERAL	Test name	GHSL	GHSL
	Analysis mode	Bireagent end point	Bireagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	280	280
	Reagent 2 (µL)	70	70
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	360 sec / cycle n° 16	360 sec / cycle n° 25	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	750	750	

xxx : User defined



GLUCOSE PAP SL

REF:
 GPSL-0490 1 x 100 mL
 GPSL-0500 6 x 100 mL
 GPSL-0507 6 x 100 mL+STD
 GPSL-0700 4 x 250 mL
 GPSL-0707 4 x 250 mL+STD

APPLICATION A15/A25	PROPOSAL
----------------------------	-----------------


For more details, see the instructions for use:
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	GPSL	GPSL
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
Reference	-	-	
Times			
Reading 1 (sec/cycle)	600 sec / cycle n° 26	600 sec / cycle n° 41	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	400	400

xxx : User defined

(09/2009)
 AA-A15/A25-GPSL-1

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C


		A15	A25
GENERAL	Test name	FEFR	FEFR
	Analysis mode	Bireagent differential	Bireagent differential
	Sample Type	serum	serum
	Units	µg/dL	µg/dL
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
Volumes	Sample (µL)	60	60
	Reagent 1 (µL)	200	200
	Reagent 2 (µL)	10	10
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
Filters	Main	560	560
	Reference	-	-
Times	Reading 1 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5
	Reading 2 (sec/cycle)	432sec / cycle n° 19	420 sec / cycle n° 29
	Reagent 2 (sec/cycle)	120 sec / cycle n° 6	120 sec / cycle n° 9
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	900	900

xxx : User defined

REF:
LACT-0100 10 x 10 mL

APPLICATION A15/A25


PROPOSAL

For more details, see the instructions for use: 
Working temperature : 37°C

		A15	A25
GENERAL	Test name	LACT	LACT
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	535	535
Reference	-	-	
Times			
Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	120	120	

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C


		A15	A25
GENERAL	Test name	LDSL	LDSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	4	4
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	24 sec / cycle n° 2	30 sec / cycle n° 3	
Reading 2 (sec/cycle)	144 sec / cycle n° 7	135 sec / cycle n° 10	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	1500	1500	

xxx : User defined

REF: LDSL-0410 2 x 62.5 mL
 LDSL-0430 4 x 62.5 mL

APPLICATION A15/A25


PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	LDSL	LDSL
	Analysis mode	Bireagent Kinetic	Bireagent Kinetic
	Sample Type	serum	serum
	Units	U/L	U/L
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
Volumes	Sample (µL)	4.5	4.5
	Reagent 1 (µL)	260	260
	Reagent 2 (µL)	65	65
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
Filters	Main	340	340
	Reference	-	-
Times	Reading 1 (sec/cycle)	72 sec / cycle n° 4	90 sec / cycle n° 7
	Reading 2 (sec/cycle)	192 sec / cycle n° 9	195 sec / cycle n° 14
	Reagent 2 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	1500	1500

xxx : User defined

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	MAGN	MAGN
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
Reference	-	-	
Times			
Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	5	5

xxx : User defined

APPLICATION A15/A25
PROPOSAL


 For more details, see the instructions for use: 

Working temperature : 37°C

		A15	A25
GENERAL	Test name Analysis mode Sample Type Units Reaction type Decimals No. of replicates Test name in patient report	PRTP Monoreagent end point serum mg/dL increasing 0 1 xxx	PRTP Monoreagent end point serum mg/dL increasing 0 1 xxx
PROCEDURE	Reading Volumes Sample (µL) Reagent 1 (µL) Reagent 2 (µL) Washing (mL) Predilution factor Postdilution factor Filters Main Reference Times Reading 1 (sec/cycle) Reading 2 (sec/cycle) Reagent 2 (sec/cycle)	Monochromatic 5 250 - 1.2 - xxx 600 - 312 sec / cycle n° 14 - -	Monochromatic 5 250 - 1.2 - xxx 600 - 300 sec / cycle n° 21 - -
CALIBRATION	Calibration type Calibrator replicates Blank replicates Calibration curve	Multiple 3 3 -	Multiple 3 3 -
	Blank absorbance Limit Kinetik Blank limit Linearity Limit	xxx - 200	xxx - 200

xxx : User defined

APPLICATION A15/A25
PROPOSAL

 For more details, see the instructions for use: 

Working temperature : 37°C

		A15	A25
GENERAL	Test name	PHOS	PHOS
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	20	20

xxx : User defined

APPLICATION A15/A25
PROPOSAL

 For more details, see the instructions for use: 

Working temperature : 37°C

		A15	A25
GENERAL	Test name	PRTB	PRTB
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	g/dL	g/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	560	560
Reference	-	-	
Times			
Reading 1 (sec/cycle)	600 sec / cycle n° 26	600 sec / cycle n° 41	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	10	10	

xxx : User defined



TRIGLYCERIDES MONO SL NEW

REF:
 TGML-0425 6 x 50 mL
 TGML-0515 6 x 100 mL
 TGML-0700 4 x 250 mL
 TGML-0427 6 x 50 mL + STD
 TGML-0517 6 x 100 mL + STD
 TGML-0707 4 x 250 mL + STD

APPLICATION A15/A25	PROPOSAL
----------------------------	-----------------

For more details, see the instructions for use:
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	TGML	TGML
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	505	505
Reference	-	-	
Times			
Reading 1 (sec/cycle)	432sec / cycle n° 19	435 sec / cycle n° 30	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	1000	1000	

xxx : User defined

(09/2009)
 AA-A15/A25-TGMLN-1



REF:	
URSL-0400	2 x 62.5 mL
URSL-0420	4 x 62.5 mL
URSL-0500	5 x 125 mL
URSL-0407	2 x 62.5 mL+ STD
URSL-0427	4 x 62.5 mL+ STD
URSL-0507	5 x 125 mL+ STD

APPLICATION A15/A25

PROPOSAL


For more details, see the instructions for use:
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	URSL	URSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	300	300
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	24 sec / cycle n° 2	30 sec / cycle n° 3	
Reading 2 (sec/cycle)	144 sec / cycle n° 7	135 sec / cycle n° 10	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	300	300	

xxx : User defined

REF:
 URSL-0400 2 x 62.5 mL
 URSL-0420 4 x 62.5 mL
 URSL-0500 5 x 125 mL
 URSL-0407 2 x 62.5 mL+ STD
 URSL-0427 4 x 62.5 mL+ STD
 URSL-0507 5 x 125 mL+ STD

APPLICATION A15/A25
PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C


		A15	A25
GENERAL	Test name	URSL	URSL
	Analysis mode	Monoreagent Kinetic	Monoreagent Kinetic
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	decreasing	decreasing
	Decimals	0	0
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	3	3
	Reagent 1 (µL)	240	240
	Reagent 2 (µL)	60	60
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	340	340
Reference	-	-	
Times			
Reading 1 (sec/cycle)	72 sec / cycle n°4	90 sec / cycle n° 7	
Reading 2 (sec/cycle)	192 sec / cycle n° 9	195 sec / cycle n° 14	
Reagent 2 (sec/cycle)	48 sec / cycle n°3	60 sec/ cycle n°5	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	300	300

xxx : User defined

REF:
 AUML-0420 6 x 50 mL
 AUML-0500 6 x 100 mL
 AUML-0700 4 x 250 mL
 AUML-0427 6 x 50 mL +STD
 AUML-0507 6 x 100 mL+STD
 AUML-0707 4 x 250 mL+STD

APPLICATION A15/A25

PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	AUML	AUML
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
Volumes	Sample (µL)	6	6
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
Filters	Main	560	560
	Reference	-	-
Times	Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21
	Reading 2 (sec/cycle)	-	-
	Reagent 2 (sec/cycle)	-	-
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
	Blank absorbance Limit	xxx	xxx
	Kinetik Blank limit	-	-
	Linearity Limit	25	25

xxx : User defined



URIC ACID SL

One reagent procedure

REF:
 AUSL-0400 2 x 62.5 mL
 AUSL-0420 4 x 62.5 mL
 AUSL-0600 5 x 125 mL

APPLICATION A15/A25	PROPOSAL
----------------------------	-----------------

For more details, see the instructions for use:
 Working temperature : 37°C


		A15	A25
GENERAL	Test name	AUSL	AUSL
	Analysis mode	Monoreagent end point	Monoreagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	No. of replicates	1	1
Test name in patient report	xxx	xxx	
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	6	6
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	-	-
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	560	560
Reference	-	-	
Times			
Reading 1 (sec/cycle)	312 sec / cycle n° 14	300 sec / cycle n° 21	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	-	-	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	25	25	

xxx : User defined

REF:
 AUSL-0400 2 x 62.5 mL
 AUSL-0420 4 x 62.5 mL
 AUSL-0600 5 x 125 mL

APPLICATION A15/A25

PROPOSAL

For more details, see the instructions for use: 
 Working temperature : 37°C

		A15	A25
GENERAL	Test name	AUSL	AUSL
	Analysis mode	Bireagent end point	Bireagent end point
	Sample Type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	No. of replicates	1	1
	Test name in patient report	xxx	xxx
PROCEDURE	Reading	Monochromatic	Monochromatic
	Volumes		
	Sample (µL)	6	6
	Reagent 1 (µL)	250	250
	Reagent 2 (µL)	62	62
	Washing (mL)	1.2	1.2
	Predilution factor	-	-
	Postdilution factor	xxx	xxx
	Filters		
	Main	560	560
Reference	-	-	
Times			
Reading 1 (sec/cycle)	360 sec / cycle n° 16	360 sec / cycle n° 25	
Reading 2 (sec/cycle)	-	-	
Reagent 2 (sec/cycle)	48 sec / cycle n° 3	60 sec / cycle n° 5	
CALIBRATION	Calibration type	Multiple	Multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
Blank absorbance Limit	xxx	xxx	
Kinetik Blank limit	-	-	
Linearity Limit	25	25	

xxx : User defined