

# MICROLAB 300

## AA-MICROLAB300-CAT-14

### A

ALBUMIN .....ALBU-4..... PRO  
ALP (DEA) SL ..... PASL-4..... PRO  
ALT / GPT 4+1 SL.....ALSL4+1-4..... PRO  
AMYLASE SL.....AMSL-6..... PRO  
AST / GOT 4+1 SL.....ASSL4+1-5..... PRO

### B

BILIRUBIN TOTAL 4+1 .....BITO-5..... PRO  
BILIRUBIN DIRECT 4+1 .....BIDI-5..... PRO

### C

CALCIUM ARSENAZO .....CALA-5..... PRO  
CHLORIDE.....CHLO-4..... PRO  
CHOLESTEROL SL .....CHSL-5..... PRO  
CHOLESTEROL HDL SL 2G.....HDLL-5..... PRO  
*NEW* HDL CHOLESTEROL .....CHDL-1..... PRO  
*NEW* CHOLESTEROL LDL SL 2G.....LDLL-5..... PRO  
LDL CHOLESTEROL.....CLDL-1..... PRO  
CHOLINESTERASE .....CHES-5..... PRO  
*UPDATE* CK - MB SL .....CMSL-5..... PRO  
CK NAC SL .....CKSL-4..... PRO  
CREATININE JAFFE .....CRCO-5..... PRO  
CREATININE PAP SL.....CRSL-5..... PRO

### G

GAMMA-GT PLUS SL.....GISL-3..... PRO  
GLUCOSE HK SL .....GHSL-5..... PRO  
GLUCOSE PAP SL.....GPSL-5..... PRO

### I

IRON FERENE .....FEFE-1..... PRO

### L

LACTATE .....LACT-4..... PRO  
LDH-L SL .....LLSL-3..... PRO  
LIPASE SL .....LPSL-2..... PRO

### M

MAGNESIUM XYLIDYL .....MAGX-1..... PRO  
MAGNESIUM XB .....MGXB-1..... PRO  
MICROPROTEIN PLUS.....PRTU-3..... PRO

### P

PHOSPHORUS .....PHOS-4..... PRO

### T

TOTAL PROTEIN PLUS.....PROB-3..... PRO  
TRIGLYCERIDES MONO SL ..TGML-4..... PRO

### U

UREA UV SL .....URSL-5..... PRO  
URIC ACID MONO SL.....AUML-5..... PRO  
URIC ACID SL.....AUSL-5..... PRO


### V

VITAMIN D.....VITD-1..... PRO

### SPECIFIC PROTEINS

ANTI STREPTOLYSINE O .....ASLO-1..... PRO  
*UPDATE* CRP IP.....ICRP-4..... PRO  
*NEW* FERRITIN .....IFRT-1..... PRO  
HAPTOGLOBIN IP .....IHAP-3..... PRO  
HBA1c .....HBAC-1..... PRO  
IgA IP .....IIGA-3..... PRO  
IgG IP .....IIGG-3..... PRO  
IgM IP .....IIGM-3..... PRO  
μALBUMIN IP .....IMAL-5..... PRO  
OROSOMUCOID IP .....IORO-3..... PRO  
PREALBUMIN IP" .....IPAL-3..... PRO  
RHEUMATHOID FACTOR .....IRFA-1..... PRO  
TRANSFERRIN IP.....ITRF-3..... PRO

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT:**

The reagent is ready for use.

**PROCEDURE:**

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).


**INSTRUMENT SETTINGS**

GENERAL	
NAME	ALBU
UNITS	g/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	620
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.000*
H-ABS-RB	0.300*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


\* This value could be optimized.

\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

**APPLICATION MICROLAB 300  
PROPOSAL**

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

**INSTRUMENT SETTINGS**

**PREPARATION OF WORKING REAGENT (WR):**



Mix 4 volumes of reagent 1 with 1 volume of reagent 2.

 **PROCEDURE:**


Pipette 1000 µL of **WR** and 20 µL of sample into a test tube.

Mix and after 30 seconds incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

GENERAL	
NAME	PASL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	405
VOLUME	500
T-FACTOR	1.000
DELAY	 30
MEAS-TIME	 30
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
FACTOR	2750


*\*\* Enter data by the operator*

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-PASL-4

## EXTENDED LINEARITY

### APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

#### PREPARATION OF WORKING REAGENT (WR):

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.




#### PROCEDURE:

Pipette 1000 µL of **WR** and 50 µL of sample into a test tube.


Mix and after 30 secondes incubation, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

#### INSTRUMENT SETTINGS

GENERAL	
NAME	ALSL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	 30
MEAS-TIME	 30
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
FACTOR	 -3333


\*\* Enter data by the operator.

 Modification from the previous version

(03/20129)

AA-MICROLAB300-ALSL4+1-4

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT:**

The reagent is ready for use.

**PROCEDURE:**


Pipette 1000 µL of reagent and 10 µL of sample into a test tube.  
Mix, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	AMSL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	405
VOLUME	500
T-FACTOR	1.000
DELAY	60
MEAS-TIME	30
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
FACTOR	9227


*\*\* Enter data by the operator*

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-AMSL-6

## EXTENDED LINEARITY

### APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

#### PREPARATION OF WORKING REAGENT (WR):

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.




#### PROCEDURE:

Pipette 1000 µL of **WR** and 50 µL of sample into a test tube.


Mix and after 30 secondes incubation, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

#### INSTRUMENT SETTINGS

GENERAL	
NAME	ASSL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	 30
MEAS-TIME	 30
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
FACTOR	 -3333

\*\* Enter data by the operator.


 Modification from the previous version

(03/2019)  
AA-MICROLAB300-ASSL4+1-5

### IMPROVED PRODUCT

REF: BITD-XXXX

#### APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

#### PREPARATION OF WORKING REAGENT:

The reagents are ready for use.

#### PROCEDURE:

##### Sample blank:

Pipette 1000 µL of reagent 1 (Direct) and 100 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

##### Sample:

Pipette 800 µL of reagent 1 (Direct) and 100 µL of sample into a test tube. Mix and wait 5 minutes at 37°C, and add 200 µL of reagent 2.

Mix and wait 50 secondes at 37°C then aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagent should be brought to the working temperature (37°C).

#### EXTRA CLEANING (after each samples)

- **1x Cleaning cycle (Bleach solution 5%)**

- **1x rinse cycle (distilled water)**

The analyzer cannot clean automatically after a sample. Please do the cleaning manually by using these 2 solutions in place of a sample blank and sample measurement.

 Modification from the previous version

#### INSTRUMENT SETTINGS

GENERAL	
NAME	BIDI
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
CALIBRATION	
METHOD	FACTOR
FACTOR	21.0


\*\* Enter data by the operator.

(03/2019)  
AA-MICROLAB300-BITD Improved-5

### IMPROVED PRODUCT

REF: BITD-XXXX

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT:**

The reagents are ready for use.

**PROCEDURE:**

Sample blank:

Pipette 1000 µL of reagent 1 (Total ) and 50 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Sample:

Pipette 800 µL of reagent 1 (Total) and 50 µL of sample into a test tube. Mix and wait 5 minutes at 37°C, and add 200 µL of reagent 2.


Mix and wait 5 minutes at 37°C then aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagents should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	BITO
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
CALIBRATION	
METHOD	FACTOR
FACTOR	31

\*\* Enter data by the operator.

 Modification from the previous version


(03/2019)  
AA-MICROLAB300-BITO-Improved Product-5



**Be careful! This application requires a filter at 660 nm.**

REF: CALA-XXXX

## APPLICATION MICROLAB 300 PROPOSAL

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT:

The reagent is ready for use.

### PROCEDURE:

Pipette 1000 µL of reagent and 20 µL of sample into a test tube.  
Mix and aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).


### INSTRUMENT SETTINGS

GENERAL	
NAME	CALA
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	660
VOLUME	500
T-FACTOR	1.000
DELAY	60
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.000*
H-ABS-RB	2.300*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-CALA-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT (WR)**

The reagent is ready for use.

**PROCEDURE:**

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).


**INSTRUMENT SETTINGS**

GENERAL	
NAME	CHLO
UNITS	mEq/L
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	505
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
L-ABS-RB	0.000*
H-ABS-RB	0.200*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-CHLO-4

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

### PREPARATION OF WORKING REAGENT

The reagent is ready for use.

### PROCEDURE:


Pipette 1000 µL of reagent and 10 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS

GENERAL	
NAME	CHSL
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	505
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
L-ABS-RB	0.000*
H-ABS-RB	0.100*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

 Modification from the previous version


\* This value could be optimized.

\*\* Enter data by the operator.

# Enter standard/calibrator value.

(03/2019)  
AA-MICROLAB300-CHSL-5

**APPLICATION MICROLAB 300  
PROPOSAL**


**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ready to use  
**Reagent 2:** Ready to use.

**PROCEDURE:**

 **Sample blank:**

Pipette 1200 µL of Reagent 1 and 10 µL of sample into a test tube.  
 Mix and after 4 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

 **Sample:**

Pipette 900 µL of Reagent 1, 10 µL of sample, mix and after 5 minutes of incubation at 37°C, add 300 µL of reagent 2 into a test tube.  
 Mix and after 4 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	HDLL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	578
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-HDLL-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ready to use  
**Reagent 2:** Ready to use.

**PROCEDURE:**

Sample blank:

Pipette 1200 µL of Reagent 1 and 12 µL of sample into a test tube.  
 Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Sample:


Pipette 900 µL of Reagent 1, 12 µL of sample, mix and after 5 minutes of incubation at 37°C, add 300 µL of reagent 2 into a test tube.  
 Mix and after 2.5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagents should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**


GENERAL	
NAME	CHDL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	578
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\*\* Enter data by the operator.  
 # Enter calibrator value.

 Modification from the previous version

(08/2020)  
 AA-MICROLAB300-CHDL-1

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ready to use

**Reagent 2:** Ready to use.

**PROCEDURE:**

 **Sample blank:**

Pipette 1200 µL of Reagent 1 and 10 µL of sample into a test tube.

Mix and after 4 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

 **Sample:**

Pipette 900 µL of Reagent 1, 10 µL of sample, mix and after 5 minutes of incubation at 37°C, add 300 µL of reagent 2 into a test tube.

Mix and after 4 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	LDLL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	578
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-LDLL-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :**  **and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ready to use

**Reagent 2:** Ready to use.

**PROCEDURE:**

Sample blank:

Pipette 1200 µL of Reagent 1 and 12 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Sample:

Pipette 900 µL of Reagent 1, 12 µL of sample, mix and after 5 minutes of incubation at 37°C, add 300 µL of reagent 2 into a test tube.

Mix and after 6.5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	CLDL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	578
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


\*\* Enter data by the operator.

# Enter calibrator value.

 Modification from the previous version

(08/2020)  
AA-MICROLAB300-CLDL-1

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT:**

**Working reagent 1 (WR1)**

Dissolve the reagent 1 with the suitable volume of distilled water as stated on the label.  
Wait 15 minutes before use.

**Working reagent 2 (WR2)**


Dissolve the reagent 2 with the suitable volume of distilled water as stated on the label.  
Wait 15 minutes before use.

**PROCEDURE:**


Pipette 1500 µL of **WR1**, 10 µL of sample and 50 µL of **WR2** into a test tube.  
Mix and aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR1** should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	CHES
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	405
VOLUME	500
T-FACTOR	1.000
DELAY	2
MEAS-TIME	90
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
CALIBRATION	
METHOD	FACTOR
FACTOR	 11730


**\*\* Enter data by the operator**

 Modification from the previous version

(08/2019)  
AA-MICROLAB300-CHES-5



## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyser manual.

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT (WR):

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.




### PROCEDURE:

Pipette 1000 µL of **WR** and 80 µL of sample into a test tube.


Mix and after 150 secondes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS


GENERAL	
NAME	CMSL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
 DELAY	30
 MEAS-TIME	300
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
 FACTOR	4286

\*\* Enter data by the operator.

 Modification from the previous version

(07/2020)  
 AA-MICROLAB300-CMSL-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :**  **and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT (WR):**      **INSTRUMENT SETTINGS**

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.

**PROCEDURE:**


Pipette 1000 µL of **WR** and 40 µL of sample into a test tube.

Mix, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

GENERAL	
NAME	CKSL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	120
MEAS-TIME	180
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
FACTOR	4127

*\*\* Enter data by the operator.*


 Modification from the previous version

(03/2019)  
AA-MICROLAB300-CKSL-4

# CREATININE JAFJE

REF: CRCO-XXXX

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

### PREPARATION OF WORKING REAGENT (WR)

Mix 1 volume of the reagent 1 with 1 volume of the reagent 2.


### PROCEDURE:


Pipette 1000 µL of **WR** and 100 µL of sample into a test tube.

Mix and aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS

GENERAL	
NAME	CRCO
UNITS	mg/dL
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	505
VOLUME	500
T-FACTOR	1.000
DELAY	10
MEAS-TIME	 30
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


 Modification from the previous version

\*\* Enter data by the operator.

# Enter standard/calibrator value.

(03/2019)  
AA-MICROLAB300-CRCO-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT:**

The reagents are ready for use.

**PROCEDURE:**

Pipette 600 µL of reagent 1 and 14 µL of sample into a test tube. Mix, wait 5 minutes at 37°C and add 200 µL of reagent 2.

Mix, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagent 1 should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	CRSL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	60
MEAS-TIME	120
NON-LIN	20
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-CRSL-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :**  **and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT (WR):**      **INSTRUMENT SETTINGS**

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.

**PROCEDURE:**

Pipette 1000 µL of **WR** and 75 µL of sample into a test tube.


Mix, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

GENERAL	
NAME	GISL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	405
VOLUME	500
T-FACTOR	1.000
DELAY	60
MEAS-TIME	30
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
CALIBRATION	
METHOD	FACTOR
FACTOR	1510*


\* This factor should be validated by using *Elical 2* calibrator.

\*\* Enter data by the operator

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-GISL-3

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :**  **and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT (WR):**

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.

**PROCEDURE:**

Pipette 1000 µL of **WR** and 10 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	GHSL
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
L-ABS-RB	0.000*
H-ABS-RB	0.850*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.


\*\* Enter data by the operator.

# Enter standard/calibrator value.

<sup>=</sup> Modification from the previous version

(03/2019)  
AA-MICROLAB300-GHSL-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT:**

The reagent is ready for use.


**PROCEDURE:**

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.  
Mix and after 10 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**


GENERAL	
NAME	GPSL
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	505
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
L-ABS-RB	0.000*
H-ABS-RB	0.250*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

 Modification from the previous version

\* This value could be optimized.  
 \*\* Enter data by the operator.  
 # Enter standard/calibrator value.

(03/2019)  
AA-MICROLAB300-GPSL-5

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT:**

**Reagent 1:** Ready to use

**Reagent 2:** Ready to use.

Before use, the reagents should be brought to the working temperature (37°C).

**CALIBRATION**

Iron Standard (REF: FEFE-0600): Ready to use

**PROCEDURE:**

Sample (Patients/Calibrator/Controls)

Pipette **480 µL of Reagent 1** and **60 µL of sample** into a test tube. Mix.

Incubate 5 minutes incubation at 37°C. Add **120 µl of Reagent 2**. Mix.

Incubate 5 minutes exactly at 37°C

Immediately after incubation, use the analyser to aspirate and start measurement.

**INSTRUMENT SETTINGS**

GENERAL	
NAME	FEFE
UNITS	µg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	578
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	0.000*
H-ABS	2.300*
REF-LOW	**
REF-HIGH	**
DECIMALS	0
L-ABS-RB	-0.100*
H-ABS-RB	2.300*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	0.3*
L-ABS	0.000*
H-ABS	2.300*
STANDARDS	
CON 1	#


\* This value could be optimized

\*\* Enter data by the operator.

# Enter standard value



## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### **PREPARATION OF WORKING REAGENT (WR):**

Dissolve the reagent 2 in the suitable volume of reagent 1.  
Wait about 15 minutes before use.


### **PROCEDURE:**

Pipette 1000 µL of **WR** and 10 µL of sample into a test tube.  
Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

### **INSTRUMENT SETTINGS**


GENERAL	
NAME	LACT
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
L-ABS-RB	0.000*
H-ABS-RB	0.300*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

 Modification from the previous version

\* This value could be optimized.  
 \*\* Enter data by the operator.  
 # Enter standard/calibrator value.

(03/2019)  
AA-MICROLAB300-LACT-4

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT (WR):      INSTRUMENT SETTINGS**

Reagents are ready to use.


**PROCEDURE:**

Pipette 1000 µL of reagent 1 (**R1**) and 35 µL of sample into a test tube. Mix, wait 3 minutes at 37°C and add 250 µL of reagent 2 (**R2**). Mix, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

GENERAL	
NAME	LLSL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	100
MEAS-TIME	180
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
FACTOR	5828

\*\* Enter data by the operator

 Modification from the previous version

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :**  **and consult the analyzer manual.**

**Working temperature :** 37°C

**PREPARATION OF WORKING REAGENT (WR):**

Reagents are ready to use.

**PROCEDURE:**


Pipette 1000 µL of reagent 1 (**R1**) and 20 µL of sample into a test tube. Mix, wait 3 minutes at 37°C and add 600 µL of reagent 2 (**R2**). Mix, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagents should be brought to the working temperature (37°C).


**INSTRUMENT SETTINGS**

GENERAL	
NAME	LPSL
UNITS	U/L
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	578
VOLUME	500
T-FACTOR	1.000
DELAY	60
MEAS-TIME	150
NON-LIN	10
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	FACTOR
FACTOR	1685

**\*\* Enter data by the operator**

 Modification from the previous version

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT (WR):**

Reagent is ready to use.

**PROCEDURE:**

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.  
Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	MAGX
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	620
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	YES
WAVELENGTH 2	505
LIMITS	
L-ABS	-0.999
H-ABS	3.000
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.300
H-ABS-RB	1.000
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-MAGX-1

**APPLICATION MICROLAB 300  
PROPOSAL**

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

**PREPARATION OF WORKING REAGENT :**

Reagent is ready to use.

**PROCEDURE:**

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.  
Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).


**INSTRUMENT SETTINGS**

GENERAL	
NAME	Mg
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	620
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	YES
WAVELENGTH 2	505
LIMITS	
L-ABS	-0.999
H-ABS	3.000
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.000*
H-ABS-RB	3.000*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	-0.999
H-ABS	3.000
STANDARDS	
CON 1	#

\* This value could be optimized.


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(04/2020)  
AA-MICROLAB300-MGXB-1

**APPLICATION MICROLAB 300  
PROPOSAL**

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

**PREPARATION OF WORKING REAGENT:**

The reagent is ready for use.


**PROCEDURE:**

Pipette 1000 µL of reagent and 20 µL of sample into a test tube.  
Mix and after 5 minutes incubation at 37°C, mix again before aspirate ,aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	PRTU
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	620
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
L-ABS-RB	0.040
H-ABS-RB	0.200
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


 Modification from the previous version

\*\* Enter data by the operator.

# Enter standard value.

(03/2019)  
AA-MICROLAB300-PRTU-3

**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT):**

The reagent is ready for use.

**PROCEDURE:**

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.  
Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).


**INSTRUMENT SETTINGS**

GENERAL	
NAME	PHOS
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.050*
H-ABS-RB	0.350*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.

\*\* Enter data by the operator.

# Enter standard/calibrator value.


 Modification from the previous version

(03/2019)  
AA-MICROLAB300-PHOS-4

# TOTAL PROTEIN PLUS

REF: PROB-XXXX

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

### PREPARATION OF WORKING REAGENT:


The reagent is ready for use.

### PROCEDURE:

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.

Mix and after 10 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).

 Modification from the previous version

<b>INSTRUMENT SETTINGS</b>	
GENERAL	
NAME	PROB
UNITS	g/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.000*
H-ABS-RB	0.100*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.

\*\* Enter data by the operator.

# Enter standard/calibrator value.

(03/2019)  
AA-MICROLAB300-PROB-3




# TRIGLYCERIDES

## MONO SL NEW

REF: TGML-XXXX

### APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

#### PREPARATION OF WORKING REAGENT:

The reagent is ready for use.

#### PROCEDURE:

Pipette 1000 µL of reagent and 10 µL of sample into a test tube.  
Mix and after 10 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).


#### INSTRUMENT SETTINGS

GENERAL	
NAME	TGML
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	505
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
L-ABS-RB	0.000*
H-ABS-RB	0.200*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-TGML-4

**APPLICATION MICROLAB 300  
PROPOSAL**

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT (WR):**

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.


**PROCEDURE:**

Pipette 1000 µL of **WR** and 10 µL of sample into a test tube.

Mix, aspirate the mixing by the analyser and start measurement.


Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	URSL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	KINETIC
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	30
MEAS-TIME	 30
NON-LIN	20
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\*\* Enter data by the operator.

# Enter standard/calibrator value.


 Modification from the previous version

(03/2019)  
AA-MICROLAB300-URSL-5

# URIC ACID MONO SL

REF: AUML-XXXX

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

### PREPARATION OF WORKING REAGENT:

The reagent is ready for use.

### PROCEDURE:

Pipette 1000 µL of reagent and 25 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagent should be brought to the working temperature (37°C).


### INSTRUMENT SETTINGS

GENERAL	
NAME	AUML
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.000*
H-ABS-RB	0.100*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-AUML-5

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

### PREPARATION OF WORKING REAGENT (WR):

Mix 4 volumes of reagent 1 with 1 volume of reagent 2.

### PROCEDURE:

Pipette 1000 µL of **WR** and 25 µL of sample into a test tube.

Mix and after 5 minutes incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the working reagent **WR** should be brought to the working temperature (37°C).


### INSTRUMENT SETTINGS

GENERAL	
NAME	AUSL
UNITS	mg/dL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
L-ABS-RB	0.000*
H-ABS-RB	0.250*
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#

\* This value could be optimized.

\*\* Enter data by the operator.


# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-AUSL-5

This application requires a filter at 700 nm.

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ready to use

**Reagent 2 :** Ready to use

### CALIBRATION:

**VITAMIN D CALIBRATOR SET :** Ready to use

### PROCEDURE:

Sample (Calibrators / Controls / Patients ):

Pipette **540 µL of Reagent 1** and **10 µL of sample** into a test tube.

Mix, wait 5 minutes at 37°C and add **135 µL of reagent 2**.

Mix, aspirate the mixing by the analyser and start measurement.

### **Note:**

**(1)** Before use, the reagents should be brought to the working temperature (37°C).


**(2)** Prime with Bleach 5% then flush with Water before each batch of Vitamin D samples and after each batch to remove latex residues.

### INSTRUMENT SETTINGS

GENERAL	
NAME	VITD
UNITS	ng/mL
R-BLANK	NO
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	TWO-POINT
WAVELENGTH	700
VOLUME	500
T-FACTOR	1.000
DELAY	24
INTERVAL	236
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	5/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#

\*\* Enter data by the operator.

# Enter standard/calibrator value.


 Modification from the previous version

(10/2019)  
AA-MICROLAB300-VITD-1

# ANTI-STREPTOLYSINE O

REF: ASLO-XXXX

**APPLICATION MICROLAB 300  
PROPOSAL**

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

**INSTRUMENT SETTINGS**

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ready to use

**Reagent 2 :** Ready to use

**CALIBRATION:**

Use **ASLO CALIBRATOR**.

**PREDILUTE CALIBRATOR, CONTROLS and PATIENTS 1:5 before use**

**PROCEDURE:**

Sample (Calibrators / Controls / Patients ):  
Pipette **600 µL of Reagent 1** and **15 µL of sample** into a test tube.

Mix, wait 5 minutes at 37°C and add **75 µL of reagent 2**.

Mix, aspirate the mixing by the analyser and start measurement.


**Note:**

Before use, the reagents should be brought to the working temperature (37°C).

GENERAL	
NAME	ASLO
UNITS	U/mL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	TWO-POINT
WAVELENGTH	546
VOLUME	500
T-FACTOR	1.000
DELAY	6
INTERVAL	130
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	1-POINT
NAME	**
FACTOR	1
REPLICATION	3
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#


\*\* Enter data by the operator.

# Enter calibrator value.

 Modification from the previous version

(10/2019)  
AA-MICROLAB300-ASLO-1

**APPLICATION MICROLAB 300  
PROPOSAL**

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ref: ICRP-6125      Ready to use  
**Reagent 2:** Ref: ICRP-5025      Ready to use

**CALIBRATION:**

**CRP IP CALIBRATOR SET**

Ref: ICRP-0043 Ready to use.  
Add a zero point (NaCl 9 g/L solution).

**PROCEDURE:**

Sample blank:

Pipette 320 µL of Reagent 1 and 20 µL of sample into a test tube.  
Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Sample:

Pipette 320 µL of Reagent 1, 32 µL of Reagent 2 and 20 µL of sample into a test tube.  
Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagents should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**


GENERAL	
NAME	ICRP
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	☞ 250
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#

\*\* Enter data by the operator.  
 # Enter standard/calibrator value.

☞ Modification from the previous version

(11/2020)  
 AA-MICROLAB300-ICRP-4

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyser manual.

Working temperature : 37°C

### INSTRUMENT SETTINGS

#### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ready to use

**Reagent 2:** Ready to use

#### CALIBRATION:

##### **FERRITIN CALIBRATOR**

Ref: IFRT-0042

Prepare a dilution range as mentioned in the Ferritin Calibrator insert sheet.

#### PROCEDURE:

Sample (Calibrators / Controls / Patients ):  
Pipette 360 µL of Reagent 1 and 40 µL of sample into a test tube.

Mix, wait 5 minutes at 37°C and add 90 µL of reagent 2.

Mix, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagents should be brought to the working temperature (37°C).

GENERAL	
NAME	FERRITIN
UNITS	ng/mL
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	TWO-POINT
WAVELENGTH	546
VOLUME	250
T-FACTOR	1.000
DELAY	6
INTERVAL	397
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	FERRITIN CAL.
NBR/ACCURACY	6/1%
REPLICATION	1
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#
CON 1	#


\*\* Enter data by the operator.

# Enter calibrator value.

(07/2020)  
AA-MICROLAB300-IFRT-1



**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyzer manual.**

**Working temperature** : 37°C

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ref: IHAP-6125 Ready to use  
**Reagent 2:** Ref: IHAP-5025 Ready to use

**CALIBRATION:**

**PROTEIN IP CALIBRATOR SET**

Ref: IPRO-0043 Ready to use.  
 Add a zero point (NaCl 9 g/L solution).

**PROCEDURE:**

Sample blank:

Pipette 840 µL of Reagent 1 and 6 µL of sample into a test tube.  
 Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Sample:

Pipette 840 µL of Reagent 1, 120 µL of Reagent 2 and 6 µL of sample into a test tube.  
 Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

**INSTRUMENT SETTINGS**

GENERAL	
NAME	IHAP
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-IHAP-3

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ready to use

**Working Reagent 2:** Pour entire contents of R2b into R2a. Mix gently (stability: 4 weeks at 2-8°C)

**Reagent 3:** Ready to use (lysing agent).

### CALIBRATION:

**HBa1C CALIBRATOR SET** (REF: HBAC-0043)

Include a zero point (NaCl, 9g/L solution).

Reconstituted Calibrators should undergo lysis pre-treatment step indicated below).

### PREPARATION (LYSIS) OF SAMPLES (PATIENTS/CALIBRATORS/ CONTROLS):

**Specimen:** Whole blood collected on EDTA. **Lysis:** Pipette **1mL of reagent R3** (lysing agent) in a plastic or glass tube. Add **20µL of a well-mixed blood sample** (calibrator, control or patient sample). Allow 5 mins for complete lysis

### PROCEDURE:

**Sample (Patients/ Calibrators/ Controls):**


Pipette **675 µL of Reagent 1** and **18 µL of lysed sample** into a test tube. Mix.

Incubate for 5 minutes at 37°C. Add **225 µL of Reagent 2**. Mix.

Incubate for 6 minutes at 37°C. Immediately after incubation, use analyser to aspirate mixture and start measurement.

### Notes:

- (1) Before use, reagents should be brought to working temperature (37°C).
- (2) Prime with Methanol then flush with Water before each batch of HBA1c samples and after each batch to remove latex residues.
- (3) If not installed by default, a 660nm filter may be purchased from ELITech and installed on instrument by a trained service person

 Modification from the previous version


### INSTRUMENT SETTINGS

GENERAL	
NAME	HbA1c
UNITS	%
R-BLANK	YES
S-BLANK	NO
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	660
VOLUME	500
T-FACTOR	1.000
DELAY	3
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	1
CALIBRATION	
METHOD	4PLL
NAME	**
NBR/ACCURACY	5/1%
REPLICATION	**
DEVIATION	0.035
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	0.0%
CON 2	#
CON 3	#
CON 4	#
CON 5	#

\*\* Enter data by the operator.

# Enter standard/calibrator value.

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ref: IIGB-6125 Ready to use

**Reagent 2:** Ref: IIGA-5025 Ready to use

### CALIBRATION:

#### **PROTEIN IP CALIBRATOR SET**

Ref: IPRO-0043 Ready to use.

Add a zero point (NaCl 9 g/L solution).

### PROCEDURE:

#### Sample blank:

Pipette 900 µL of Reagent 1 and 5 µL of sample into a test tube.

Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

#### Sample:

Pipette 900 µL of Reagent 1, 150 µL of Reagent 2 and 5 µL of sample into a test tube.

Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS

GENERAL	
NAME	IIGA
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-IIGA-3

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### **PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ref: IIGB-6125      Ready to use  
**Reagent 2:** Ref: IIGG-5025      Ready to use

### **CALIBRATION:**

#### **PROTEIN IP CALIBRATOR SET**

Ref: IPRO-0043 Ready to use.  
 Add a zero point (NaCl 9 g/L solution).

### **PREPARATION OF SAMPLES AND CONTROLS**

Dilute samples and controls at 1/10 in NaCl 9g/L solution.

### **PROCEDURE:**


#### Sample blank:

Pipette 500 µL of Reagent 1 and 10 µL of diluted sample into a test tube.  
 Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

#### Sample:

Pipette 500 µL of Reagent 1, 50 µL of Reagent 2 and 10 µL of diluted sample into a test tube.  
 Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagents should be brought to the working temperature (37°C).

 Modification from the previous version


### **INSTRUMENT SETTINGS**

GENERAL	
NAME	IIGG
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	250
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#

\*\* Enter data by the operator.

# Enter standard/calibrator value.

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ref: IIGB-6125      Ready to use  
**Reagent 2:** Ref: IIGM-5025      Ready to use

### CALIBRATION:

#### **PROTEIN IP CALIBRATOR SET**

Ref: IPRO-0043 Ready to use.  
 Add a zero point (NaCl 9 g/L solution).

### PROCEDURE:

#### Sample blank:

Pipette 900 µL of Reagent 1 and 5 µL of sample into a test tube.  
 Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

#### Sample:

Pipette 900 µL of Reagent 1, 130 µL of Reagent 2 and 5 µL of sample into a test tube.  
 Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS

GENERAL	
NAME	IIGM
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
 AA-MICROLAB300-IIGM-3

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ref: IMAL-6125 Ready to use  
**Reagent 2:** Ref: IMAL-5025 Ready to use

### CALIBRATION:

#### μALBUMIN IP CALIBRATOR SET

Ref: IMAL-0043 Ready to use.  
Add a zero point (NaCl 9 g/L solution).

### PROCEDURE:

#### Sample blank:

Pipette 300 μL of Reagent 1 and 20 μL of sample into a test tube.  
Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

#### Sample:

Pipette 300 μL of Reagent 1, 50 μL of Reagent 2 and 20 μL of sample into a test tube.  
Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS

GENERAL	
NAME	IMAL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	250
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	2
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-IMAL-5

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ref: IORO-6125 Ready to use

**Reagent 2:** Ref: IORO-5025 Ready to use

### CALIBRATION:

#### PROTEIN IP CALIBRATOR SET

Ref: IPRO-0043 Ready to use.

Add a zero point (NaCl 9 g/L solution).

### PROCEDURE:

#### Sample blank:

Pipette 750 µL of Reagent 1 and 5 µL of sample into a test tube.

Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

#### Sample:

Pipette 750 µL of Reagent 1, 100 µL of Reagent 2 and 5 µL of sample into a test tube.

Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS

GENERAL	
NAME	IORO
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#


\*\* Enter data by the operator.

# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-IORO-3

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

**Working temperature** : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ref: IPAL-6125 Ready to use

**Reagent 2:** Ref: IPAL-5025 Ready to use

### CALIBRATION:

#### **PROTEIN IP CALIBRATOR SET**

Ref: IPRO-0043 Ready to use.

Add a zero point (NaCl 9 g/L solution).

### PROCEDURE:

#### Sample blank:

Pipette 550 µL of Reagent 1 and 5 µL of sample into a test tube.

Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

#### Sample:

Pipette 550 µL of Reagent 1, 50 µL of Reagent 2 and 5 µL of sample into a test tube.

Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

### INSTRUMENT SETTINGS

GENERAL	
NAME	IPAL
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	250
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#

\*\* Enter data by the operator.


# Enter standard/calibrator value.

 Modification from the previous version

(03/2019)  
AA-MICROLAB300-IPAL-3



**APPLICATION MICROLAB 300  
PROPOSAL**

**For more details, see the Instructions for use :  and consult the analyser manual.**

**Working temperature** : 37°C

**INSTRUMENT SETTINGS**

**PREPARATION OF WORKING REAGENT :**

**Reagent 1:** Ready to use

**Reagent 2:** Ready to use

**CALIBRATION:**

**RF CALIBRATOR SET**

Ref: IRFA-0042 Ready to use.

Add a zero point (NaCl 9 g/L solution).

Prepare a dilution range as mentioned in the RF Calibrator insert sheet.

**PROCEDURE:**

Sample blank:

Pipette 400 µL of Reagent 1 and 5 µL of sample into a test tube.

Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Sample:

Pipette 400 µL of Reagent 1, 100 µL of Reagent 2 and 5 µL of sample into a test tube.

Mix and after 2 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.


Note: Before use, the reagents should be brought to the working temperature (37°C).

GENERAL	
NAME	RF
UNITS	U/mL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	620
VOLUME	250
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	0.000
H-ABS	3.000
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	RF CALIBRATOR
NR/ACCURACY	6/1%
REPLICATION	1
DEVIATION	-
L-ABS	0.000
H-ABS	3.000
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#

\*\* Enter data by the operator.

# Enter calibrator value.

## APPLICATION MICROLAB 300 PROPOSAL

For more details, see the Instructions for use :  and consult the analyzer manual.

Working temperature : 37°C

### PREPARATION OF WORKING REAGENT :

**Reagent 1:** Ref: ITRF-6125 Ready to use

**Reagent 2:** Ref: ITRF-5025 Ready to use

### CALIBRATION:

#### **PROTEIN IP CALIBRATOR SET**

Ref: IPRO-0043 Ready to use.

Add a zero point (NaCl 9 g/L solution).

### PREPARATION OF SAMPLES AND CONTROLS

Dilute samples and controls at 1/20 in NaCl 9g/L solution.

### PROCEDURE:

#### Sample blank:

Pipette 880 µL of Reagent 1 and 12 µL of diluted sample into a test tube.


Mix and after 1 minute incubation at 37°C, aspirate the mixing by the analyser and start measurement.

#### Sample:

Pipette 880 µL of Reagent 1, 40 µL of Reagent 2 and 12 µL of diluted sample into a test tube.

Mix and after 5 minutes of incubation at 37°C, aspirate the mixing by the analyser and start measurement.

Note: Before use, the reagents should be brought to the working temperature (37°C).

 Modification from the previous version

### INSTRUMENT SETTINGS

GENERAL	
NAME	ITRF
UNITS	mg/dL
R-BLANK	NO
S-BLANK	YES
REPLICATION	1
SETTINGS	
MODE	END-POINT
WAVELENGTH	340
VOLUME	500
T-FACTOR	1.000
DELAY	2
BI-CHROM	NO
WAVELENGTH 2	-
LIMITS	
L-ABS	**
H-ABS	**
REF-LOW	**
REF-HIGH	**
DECIMALS	0
CALIBRATION	
METHOD	MULTI-POINT
NAME	**
NBR/ACCURACY	6/1%
REPLICATION	**
DEVIATION	**
L-ABS	**
H-ABS	**
STANDARDS	
CON 1	#
CON 2	#
CON 3	#
CON 4	#
CON 5	#
CON 6	#

\*\* Enter data by the operator.

# Enter standard/calibrator value.

(03/2019)  
AA-MICROLAB300-ITRF-3